



0.6x - 5x Zoom Stereo 1 / 8.3 Stereo Optical Microscope A23.0404

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC
- Certification: CE, Rohs
- Model Number: A23.0404
- Minimum Order Quantity: 1
- Price: FOB \$100~1000, Depend on Order Quantity
- Packaging Details: Carton Packing
- Delivery Time: 5-30
- Payment Terms: L/C, T/T, Western Union
- Supply Ability: 10000



Product Specification

- Highlight: stereo zoom microscopes,
stereo dissecting microscope



Product Description

0.6x - 5x Zoom Stereo 1 / 8.3 Stereo Optical Microscope A23.0404

A23.0404 Zoom Stereo Microscope		B	B	T	T
		I			
Head	Binocular head, 45° Inclined interpupillary distance: 48-75mm, diopter adjustment: ±5 diopter	o	o		
	Trinocular head, 45° Inclined interpupillary distance: 48-75mm, diopter adjustment: ±5 diopter			o	o
Eyepiece	SWF10x/22mm				
Zoom Ratio	1:8.3				
Objective	0.6x-5x				
Zoom Range					
Stand	Track stand				
Working Distance	95mm				
Illumination	No light	o	o		
	Top halogen lamp 6v/15w, adjustable brightness	o	o		
	Bottom fluorescent lamp 5w	o	o		
Voltage Input	110V/60HZ or 220V/50HZ				
A23.0404 Zoom Stereo Microscope Optional Accessories					
Eyepiece	5x	A51.04			
		01-05			
	10x	A51.04			
		01-10			
	15x	A51.04			
		01-15			
	20x	A51.04			
Auxiliary Objective		01-20			
	25x	A51.04			
		01-25			
	30x	A51.04			
		01-30			
	0.5x	A52.04			
		05-05			
Photo Adapter	6.3x	A55.04			
		05-63			
	2.5x	A55.04			
CCD Adapter		05-25			
	0.5x	A55.04			
		13-05			
Base	1x	A55.04			
		13-10			
	Pole stand, no light	A55.04			
		02			
CCD	Pole stand, reflected and transmitted illumination	A55.04			
		03			
	Pole stand, reflected and transmitted illumination, LED light (110V-240V)	A55.04			
Eye Adapter		04			
		A55.04			
Polarizing Attachment		16			
	Polarizer	A57.04			
		0			



A24.0404-B

A24.0404-BI



A24.0404-T

A24.0404-TI



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China