



10x - 40x Trinocular Stereo Optical Microscope Track Stand and Pole Stand A22.0301

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC
- Certification: CE, Rohs
- Model Number: A22.0301
- Minimum Order Quantity: 1
- Price: FOB \$100~1000, Depend on Order Quantity
- Packaging Details: Carton Packing
- Delivery Time: 5-30
- Payment Terms: L/C, T/T, Western Union
- Supply Ability: 10000



Product Specification

- Highlight: stereo zoom microscopes,
stereo binocular microscopes



More Images



for more products please visit us on cnoec.com

Product Description

10x - 40x Trinocular Stereo Optical Microscope Track Stand and Pole Stand A22.0301

A22.0301 stereo microscope, suitable for observe and research in the field of teaching biological anatomy, electronic engineering, geology & mine, farming and forestry of plant protection, penal reconnaissance, jewelry appraisal and so on.

--Objective 1x/2x, 1x/3x, 2x/4x optional

--SWF10x/2mm eyepiece optional

--Different stand with or without lightsource for choice

Model	44	54	47	57	46	56	49	59
Head	Binocular	Trinocular	Binocular	Trinocular	Binocular	Trinocular	Binocular	Trinocular
Stand	Pole Stand		Track stand		Pole Stand		Track stand	
	D, E, F Stand Optional							
Light	Up/Bottom light source				No light source			
		Specification						
Head		Binocular / Trinocular, inclined 45°, interpupillary distance 55-75mm						
Eyepiece		WF10x//20mm						
Objective		1x/2x, 1x/3x, 2x/4x						
Working Distance		1x/2x 110mm, 1x/3x 102mm, 2x/4x 90mm						
Viewing Field		23mm-5.7mm						
Working Stage		Round Working Stage Φ95mm, filter Φ40mm						
Illumination		Up light 12V 10W halogen lamp, brightness adjustable						
		Bottom light fluorescent lamp, brightness adjustable						
Optional Accessories							Item No.	
Eyepiece			SWF10x/23mm				A51.0306	
			WF5x				A51.0305-05	
			WF12.5x				A51.0305-12.5	
			WF15x				A51.0305-15	
			WF20x				A51.0305-20	
			WF25x				A51.0305-25	
Auxiliary Objective			0.5x				A52.0380-0.5	
			0.75x				A52.0380-0.75	
			1.5x				A52.0380-1.5	
			2.0x				A52.0380-2	
Adapter			CCD adapter				A55.0302	
			Photo adapter				A55.0301	
			Digital camera adapter				A55.0303	
Condenser			Dark field condenser				A57.0310	

A22.0301-44
Binocular, Pole Stand

www.cnoec.com



www.cnoec.com




A22.0301-54
Trinocular, Pole Stand

www.cnoec.com



A22.0301-47
Binocular, Track Stand

A22.0301-57
Trinocular, Track Stand



D Base



E Base



F Base



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China

