



## Vertical Monocular Head Biological Microscope with 5 holes diagram Condenser A11.0306

### Our Product Introduction

#### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.0306
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



#### Product Specification

- Applications: Student
- Certification: CE
- Drawtube: Monocular
- Illumination: LED
- Product Name: Microscope
- Theory: Biological Microscope
- Highlight: **high powered microscope ,  
biological compound microscope**



## Product Description

### Vertical Monocular Head Biological Microscope with 5 holes diagram Condenser A11.0306

**A11.0306 Student Microscope**, an ideal tool or gift scope for children and student, with full accessories kit all in one..




--40x~1600x Magnification

--Full Metal Body for Steady Observation

--Compact All in One Design for Easy Moving

Specific ation	A11.0306 Biological Microscope	AB
Head	Vertical Monocular Head	
Eyepiece	H10x	
	H16x	
Nosepiece	Triple	●
	Quadruple	●
Objective	Achromatic Objective 4x/NA 0.1	
	Achromatic Objective 10x/NA 0.25	
	Achromatic Objective 40x(s)/NA 0.65	
	Achromatic Objective 100x(s,oil)/NA.1.25	●
Stage	Fixed stage size 115x125mm,stage with clip	●
	Fixed stage size 115x125mm, moving scale,moving range:60x30mm, minimum vernier division 0.1mm	●
Condenser	N.A.0.65 with 5 holes diagram	●
	N.A.1.25 Abbe condenser, up/down range 15mm	●
Focusing	Coaxial coarse & fine focusing knob, coarse focusing range 35mm, fine focusing division 0.002mm	
Filter	Blue filter , Yellow filter, Ground glass	
Light Source	Plano-Concave mirror: Φ 50mm	



 0086 13911110627    [sale@optoedu.com](mailto:sale@optoedu.com)    [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China