



0.65x-4.5x Zoom Stereo 1 / 7 Stereo Optical Microscope Tinocular and Binocular A23.0403

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC
- Certification: CE, Rohs
- Model Number: A23.0403
- Minimum Order Quantity: 1
- Price: FOB \$100~1000, Depend on Order Quantity
- Packaging Details: Carton Packing
- Delivery Time: 5-30
- Payment Terms: L/C, T/T, Western Union
- Supply Ability: 10000



Product Specification

- Highlight: stereo zoom microscopes, stereo dissecting microscope



Product Description

Tinocular and Binocular 0.65x-4.5x Zoom Stereo 1 / 7 Stereo Optical Microscope A23.0403

This microscope is a popular design, as well as excellent optical performance

A23.0403 Zoom Stereo Microscope		B B I T T I	
Head	Binocular head, 45°/60° Inclined interpupillary distance: 52-75mm, diopter adjustment: ±5 diopter	o	o
	Trinocular head, 45°/60° Inclined interpupillary distance: 52-75mm, diopter adjustment: ±5 diopter		o o
Eyepiece	SWF10x/22mm		
Zoom Ratio	1:7		
Objective	0.65x-4.5x		
Zoom Range			
Stand	Track stand		
Working Distance	103mm		
Illumination	No light	o	o
	Top halogen lamp 6v/15w, adjustable brightness	o	o
	Bottom fluorescent lamp 5w	o	o
Voltage Input	110V/60HZ or 220V/50HZ		
A23.0403 Zoom Stereo Microscope Optional Accessories			
Eyepiece	5x	A51.040	1-05
	15x	A51.040	1-15
	20x	A51.040	1-20
	25x	A51.040	1-25
	30x	A51.040	1-30
		A51.040	1-30
Auxiliary Objective	0.5x	A52.040	6-05
	0.75x	A52.040	6-75
	1.5x	A52.040	3-15
	2x	A52.040	3-20
Photo Adapter	6.3x	A55.040	5-63
	2.5x	A55.040	5-25
CCD Adapter	0.5x	A55.041	3-05
	1x	A55.041	3-10
Polarizing Attachment	Polarizer	A5P.040	1-1



A23.0403-B
Binocular Zoom Stereo Microscope





Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China