



Dia 95mm Microscope Accessores Digital Working StageA54.0401

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A54.0401
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Microscope
- Certification: CE|Rohs
- Product Name: Microscope Accessories
- Resolution: 2.0M
- Senser: CMOS
- Pixel Size: 2.75*2.75um
- Highlight: **compound microscope parts, parts of the microscope**

for more products please visit us on cnoec.com

Product Description

Dia 95mm Digital Working Stage Microscope Accessory e A54.0401

This stage is suitable to fix onto stereo microscope stand, which has 95mm hole for plate. it is flexible to move and operate.

A54.0401 Digital Working Stage	
Specification of Stage:	--Stage size:140x140x65mm (LxWxH) --Net Weight:2.5kgs --Connetion size for stand: 95mm
Specification of Micrometer	--Measuring:0-25mm/0-1" --Resolution:0.001mm/0.00005" --Accuracy:0.002mm --Power supply:one silver-oxide button battery SR44(1.5V) --Operating temperature:0-+40° --Storage temperature:-20-70° --No influence within relatvie humidity<80%



A54.0401 Digital Working Stage



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627 ✉ sale@optoedu.com 🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China