



USB3.0 16MP Digital Microscope Cameras Drive Mode With 16 Bits A59.3515

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A59.3515
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Microscope
- Certification: CE|Rohs
- Product Name: Microscope Accessories
- Resolution: 16.0MP
- Output: USB 3.0
- Trigger: External Trigger/software Trigger
- Highlight: **microscope accessory, parts of the microscope**



More Images



Product Description

USB3.0 16MP Digital Microscope Cameras Drive Mode With 16 Bits A59.3515

- Max resolution can reach 16MP
- The output interface is USB 3.0
- Standard C-mount

A59.3515 Digital Camera

Senser	COMS
Resolution	4632*3488@8,1080P@60,720P@60
Pixel Size	1.335um*1.335um
Data Output	16 Bits
Data Format	YUV422
Partial Image	ROI
Image Processing	GAMMA,LUT, Aaturation, Acutance
Eletronic Shutter(us)	20-100000(full resolution),step:1
Gain(db)	0-36,step:0.1
White Balance(db)	0-12,step:0.1
Gamma	0.5-4.00
Digial Interface	USB 3.0
GPIO	8pin, HIROSE HR25:PWR,1XOPI,1XOPO,BIOX2
Trigger	External trigger/software trigger
Power	DC5V,<4.5W, GPIO/USB
Len Interface	C-mount
Camera Control	Drive Mode
Temperature	0-50




A59.3515 Digital Camera, USB 3.0, 16MP





The Lens Is Optional



www.cnoec.com

A59.3515 Digital Camera, USB 3.0, 16MP



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China