



40x - 400x Laboratory Research Polarization Microscope Binocular Halogen Lamp A15.2604

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A15.2604
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Application: Research|Laboratory
- Certification: CE|Rohs
- Drawtube: Binocular
- Illumination: Halogen Lamp
- Magnification: 40X - 600X
- Theory: Polarizing Light Microscope
- Highlight: **highest magnification microscope, compound light microscope**



More Images





OPTO-EDU (BEIJING) CO., LTD.

F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China
Tel: +8610 88696020 Fax: +8610 88696085

A15.2604

Polarizing Microscope





Monocular Head, WF 10X/18mm, With Cross Micrometer 0.1mm

Head & Eyepiece

Monocular Head, Inclined 30°, Rotatable 360°
WF 10X/18mm, With Cross Micrometer 0.1mm



Analyzer Unit, 0-90° , It Can Be Moved Out Of
The Optical Path For Single Polarizing Observe

Analyzer Unit

0-90°, It Can Be Moved Out Of The Optical Path For Single Polarizing Observe


Bertrand Lens

It Can Be Moved Out Of The Optical Path



Strain-free Achromatic Objective 4x,10x 40x(S)


Quadruple Nosepiece



Round Stage, Diameter 142mm, Rotatable 360°

Working Stage

Round Stage, Diameter 142mm, Rotatable 360°, Division 1°, Vernier division 6, Lockable



Abbe N.A. 1.25 Condenser



12V/20W Halogen Lamp

Condenser & Illumination

Abbe N.A. 1.25, 12V/20W Halogen Lamp, With Blue Filter On Light Collector



Coaxial Coarse & Fine Focusing



Powered By Safety Low Voltage Charger, Input Wide Voltage 100V-240V, Output 5V1A, Support Power Bank Supply For Outdoor Use

5W LED Light Source Powered By Safety Low Voltage Charger

Input Wide Voltage 100V-240V, Output 5V1A, Support Power Bank Supply For Outdoor Use

Items	A15.2604 Polarizing Microscope
Head	Monocular Head, Inclined 30°, Rotatable 360°
Eyepiece	WF 10X/18mm, With Cross Micrometer 0.1mm

Objective	Strain-free Achromatic Objective 4x,10x 40x(S)	
Nosepiece	Quadruple	
Analyzer Unit	0-90°, It Can Be Moved Out Of The Optical Path For Single Polarizing Observe	
Bertrand Lens	It Can Be Moved Out Of The Optical Path	
Polarizing	Polarizer Unit Built-in Condenser	
	Analyzer Unit Insert to Polarizing Media, Rotatable 0-90°, Can Be Moved Out Of Optical Path	
	Polarizing Media, Between Head & Body, Including Compensator:	
	--Bertrand Lens, Can Be Moved Out Of Optical Path	
	--λ Slip (First Class Red)	
	--λ/4 Slip	
Stage	Round, Diameter 142mm, Rotatable 360°, Division 1°, Vernier division 6, Lockable	
Condenser	Abbe N.A. 1.25	
Polarizing Unit	Can Be Moved Out Of The Optical Path	
Illumination	12V/20W Halogen Lamp, With Blue Filter On Light Collector	
Optional Accessory		
Eyepiece	WF16x/13mm	A51.2602-1613
Objective	Strain-free Achromatic Objective 20x	A5P.2629-20
	Strain-free Achromatic Objective 60x	A5P.2629-60
Compensator	Quartz Wedge I~IV Class	A5P.2615-17
Stage	Attachable Polarizing Mechanical Stage, Moving Range 40x40mm	A5P.2616
Illumination	3W LED Illumination, Powered By Safety Low Voltage Charger, Input Wide Voltage 100V-240V, Output 5V1A, Support Power Bank Supply For Outdoor Use	A56.2613-3W
	3W LED Illumination (Rechargeable), Powered By Safety Low Voltage Charger, Input Wide Voltage 100V-240V, Output 5V1A, Support Power Bank Supply For Outdoor Use	A56.2651-3WR

Product Accessories



A5P.2616 Attachable Polarizing Mechanical Stage, Moving Range 40x40mm



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China