



Binocular / Trinocular Stereo Optical Microscope 0.7 - 4.5x Zoom Stereo 1:6.5 A23.0906-S1

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC
- Certification: CE, Rohs
- Model Number: A23.0906-S1
- Minimum Order Quantity: 1
- Price: FOB \$100~1000, Depend on Order Quantity
- Packaging Details: Carton Packing
- Delivery Time: 5-30
- Payment Terms: L/C, T/T, Western Union
- Supply Ability: 10000



Product Specification

- Highlight: stereo zoom microscopes,
stereo dissecting microscope

Product Description

Binocular / Trinocular Stereo Optical Microscope 0.7 - 4.5x Zoom Stereo 1:6.5 A23.0906-S1

- 0.7~4.5x Zoom Lens, with WF10x/20mm Eyepiece,
- Total Optical Magnification 7~45x
- Single Arm Universal Pole Stand

A23.0906-S1 Zoom Stereo Microscope, 0.7x-4.5x Specification

		A23.0906-S1	A23.0906-S1T
Head	Binocular, 45° Inclined, 360° Rotatable, +/-5° Diopter Adjustable, Interpupillary Adjustment 54~75mm	o	
	Trinocular, 45° Inclined, 360° Rotatable, +/-5° Diopter Adjustable, Interpupillary Adjustment 54~75mm		o
Zoom Range	0.7~4.5x		
Zoom Ratio	1:6.5		
Eyepiece	WF10x/20mm, Both with Eye-cup		
Magnification	7~45x, Up to 3.5~180x with Auxiliary Lens		
Working Distance	100mm, With Standard 1x Lens		
Field of View	1~95.2mm		
Stand	Single Arm Universal Pole Stand STL1, Height 308mm, Focusing Pole Height 104mm, Level Move Range 235mm, Base 230*230*40mm,		
Illumination	No Light Source		

A23.0906-S1 Zoom Stereo Microscope, 0.7x-4.5x Optional Accessories

Eyepiece	WF15x/15mm	A51.0922-15
	WF20x/10mm	A51.0922-20
	0.5x	A52.0927-05
	0.75x	A52.0927-75
Auxiliary Lens	1.5x	A52.0927-15
	2.0x	A52.0927-20



A23.0906-S1 Zoom Stereo Microscope, 0.7~4.5x, With Universal Stand



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China