



A43.1911 Video Colposcope Digital Continuous Zoom 3x-36x

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE,
- Model Number: A43.1911
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- FHD Camera: Sony 1/2.8" Exmor R CMOS
- Max. Resolution: 1920*1080P
- Objective Focal Length: F=300mm (M58x0.75mm)
- Total Magnification: 3X-36X Zoom
- Working Distance: 200mm-350mm
- Focusing Mode: Auto And Manual
- Video Output Interface: HDMI
- Video System: Monitor 24" LCD, Monitor Holder, HDMI-USB3.0 Type C Video Captuer Box
- Floor Stand: Universal Stand,
- Microscope Height: 76mm-130mm
- Highlight: **Continuous Zoom Video Colposcope, 36x Video Colposcope, Continuous Zoom Digital Video Colposcope**

Product Description

A43.1911 Video Colposcope Digital Continuous Zoom 3x-36x

Table Of Specification For A43.1911 & A43.1912 Series Colposcope		
Item	A43.1911	A43.1912
FHD Camera	Sony 1/2.8" Exmor R CMOS	
Max. Resolution	1920*1080P	
Objective Focal Length	F=300mm (M58x0.75mm)	
Total Magnification	3X-36X Zoom	
Working Distance	200mm-350mm	
Focusing Mode	Auto And Manual	
Video Output Interface	HDMI	
Menu	Focusing, White Balance, Exposure, Color, ICR etc.	
Illumination System	Coaxial Illumination, 10W LED, Brightness Adjustable, Illumination>10000lx	
Power Supply	Power Supply Input: AC110~240V, 50/60Hz, With Power Supply Adapter Output DC12V3A, 36W,	
Filter	Built-in Green Filters	
Video System	Monitor 24" LCD, Monitor Holder, HDMI-USB3.0 Type C Video Captuer Box	
Floor Stand	Universal Stand,	Straight Vertical Stand,
	Five-star Base With Casters	Height Adjustable Range 24cm, Swing Angle $\pm 8^\circ$
Microscope Height	76mm-130mm	86mm-110mm
Optional Accessories	Vaginal Microscope Image Analysis & Reporting Software, Monitor, Objective (F=230mm & F=280mm)	



OPTO-EDU (BEIJING) CO., LTD.

OPTO-EDU


F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China
Tel:+8610 88696020 Fax:+8610 88696085

A43.1911

Digital Video Colposcope, Continous Zoom 3x-36x, Universal Stand



Digital Camera Sony 1/2.8" Exmor R CMOS, Resolution 1920*1080P, Objective F=300mm, Total Magnification 3x-36x Zoom, W.D.200-350mm, Auto/Manual Focus, HDMI Output







A43.1911



A43.1912



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China