



## USB Digital Hand Held Digital Microscope

### Our Product Introduction

#### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A34.4161
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



#### Product Specification

- Image Sensor: CMOS Sensor
- Model Number: A34.4161
- Application: Research
- Certification: CE
- Magnification: 800x
- Theory: Handheld Digital Microscope
- Output: USB
- Highlight: **Portable Digital Microscope,  
Hand held Digital Microscope**

for more products please visit us on [cnoec.com](http://cnoec.com)

## Product Description

800 X 2.0M USB Digital Hand Held Digital Microscope Inspection A34.4161

### Description:

New Product in 2015!

It can be easily used in different fields including:

- 1.Textile industrial for textile inspection
- 2.Printing inspection
- 3.Industrial inspection:PCB,Priction machinery
- 4.Educational Purpose
- 5.Hair Examination
- 6.Skin Examination
- 7.Microbiological Obervation
- 8.Jewelry&Coin(Colletions)inspection
- 9.Visual Assistance

### A34.4161 Specification

Color of Body	Black
Magnification	20 X ~ 800X ( Manual )
Focus Range	0mm ~ 40mm
Still Image Resolution	1600*1200, 640*480, 480*320, 320*240
Lens	Microscope Lens , M12
Resolution	2.0M ( Default 0.3M )
Frame Rate	30 f/s Under 600 Lux Brightness
Digital Zoom	5X Sequence Mode
Image Sensor	Color CMOS Sensor
Controller	High Speed DSP( Driver Free available )
Power Supply	USB Port ( 5V DC )
Interface	USB 2.0
Illumination	White Light 8pcs LED Manual adjustment
Snap Shot	Software and Hardware
Measure Tool	Micro-Measure(It can measure Length and diameters, angle, the perimeter, area, etc) (Language:English,Chinese Simplified)
Driver Language	English,Chinese
Operation System	Windows XP ,Vista,Windows 7, 32 bit and 64 bit
System Requirement	Pentium Computer with 700M Hz & Above , 20M HD Space CD ROM Driver, 64MB RAM , Direct X VGA Card
Accesseries	1.5M USB cable, CD disk,Metal Fix Stand, color box, calibration ruler, user manual



Magnification 20x-800x



U Fix Stand



USB interface



8pcs LED



Circuit Board



Money



Texture



Skin



Printing



**Opto-Edu (Beijing) Co., Ltd.**



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China