



Tinocular / Binocular 1 / 6.5 Stereo Optical Microscope 0.7 - 4.5x Zoom Stereo A23.5601

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC
- Certification: CE, Rohs
- Model Number: A23.5601
- Minimum Order Quantity: 1
- Price: FOB \$100~1000, Depend on Order Quantity
- Packaging Details: Carton Packing
- Delivery Time: 5-30
- Payment Terms: L/C, T/T, Western Union
- Supply Ability: 10000



Product Specification

- Highlight: stereo dissecting microscope,
stereo binocular microscopes



Product Description

Tinocular / Binocular 1 / 6.5 Stereo Optical Microscope 0.7 - 4.5x Zoom Stereo A23.5601

- 0.7~4.5x Zoom Lens, with WF10x/20mm Eyepiece,
- Total Optical Magnification 7~45x
- Single Arm Universal Pole Stand

A23.5601 Zoom Stereo Microscope Specification				B	T		
Head	Binocular Head, 45° inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 55-75mm			o			
	Trinocular Head, 45° inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 55-75mm				o		
Eyepiece	WF10x/20mm, High Eyepoint						
Zoom Lens	0.7~4.5x						
Zoom Ratio	1:6.5						
Magnification	7x-45x,						
Working Distance	100mm						
Illumination	Double Pipe 3W LED Light Source, Brightness Adjustable						
Stand	Pillar Stand						
Power	110V/220V						
A23.5601 Zoom Stereo Microscope Optical Data Sheet							
Objective	Standard 1.0x	Auxiliary 0.5X		Auxiliary 1.5X		Auxiliary 2X	
W.D.	100mm	165mm		45mm		30mm	
Eyepiece	Mag	View field	Mag	View field	Mag	View field	
10X/20mm	7.0X	28.6mm	3.5X	57.1mm	10.5X	19mm	14X
	45.0X	4.4mm	22.5X	8.9mm	67.5X	3mm	90X
15X/15mm	10.5X	21.1mm	5.25X	42.8mm	15.75X	14.3mm	21X
	67.5X	3.3mm	33.75X	6.7mm	101.25X	2.2mm	135X
20X/10mm	14.0X	14.3mm	7.0X	28.6mm	21.0X	9.5mm	28X
	90.0X	2.2mm	45.0X	4.4mm	135.0X	1.5mm	180X
							1.1mm
A235601 Zoom Stereo Microscope Optional Accessories					Cata. No.		
Head	Binocular Head + WF10x/20mm						
	Trinocular Head + WF10x/20mm						
Eyepiece	WF10x/20mm, High-eyepoint, Outside Dia. 30mm						
	WF15x/15mm, High-eyepoint, Outside Dia. 30mm						
	WF20x/10mm, High-eyepoint, Outside Dia. 30mm						
	WF10x/20mm,with Reticle						
	WF15x/15mm,with Reticle						
Auxiliary Lens	WF20x/10mm,with Reticle						
	0.5X,W.D. 165mm						
	0.75X,W.D. 105mm						
	1.5X,W.D. 45mm						
CCD Adapter	2.0X,W.D. 30mm						
	0.50x, For 1/2" CCD						
	0.35x, For 1/3" CCD						





Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China