



## Compensation Binocular Phase Contrast Microscope with Kohler 6V 20W Illumination A19.0209

### Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

#### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE,
- Model Number: A19.0209
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month





#### Product Specification

- Certification: CE
- Drawtube: Binocular
- Illumination: 6V 20W
- Theory: Phase Contrast Microscope
- Condensor: Kohler N.A. 1.25
- Nosepiece: Backward Quadruple
- Highlight: **phase contrast light microscope, inverted phase contrast microscope**



#### More Images



## Product Description

### Compensation Binocular Phase Contrast Microscope with Kohler 6V 20W Illumination A19.0209


- Ergonomic Design, Easy Operation, Outstanding Optical Performance, Clear and Bright Image
- Specially Developed for Phase Contrast Observation, High Contrast, Highlight the Cell Outline and Internal Structure
- High Precision Focusing Device and Illumination Adjustment


A19.0209 Phase Contrast Microscope	
Head	Compensation binocular head inclination 30°;
Eyepiece	WF10X(Φ18mm)
Nosepiece	Backward Quadruple Nosepiece, Ball Bearing Inner Locating
Phase Contrast Set	Phase contrast Achromatic Objective:
	10x/0.25 PHP
	20x/0.40 PHP(Spring)
	40x/0.65 PHP(Spring)
	100x/1.25 PHP (Spring,oil)
	Centering telescope
	Turnplate phase contrast condenser
Working Stage	Double layer mechanical, size 135x125mm, moving range 75x35mm
Focusing	Coaxial coarse/fine focus system, Minimum division of fine focusing: 2um, With tensional adjustable and up stop, Base bearing black plastic cap
Condenser	Abbe N.A.1.25 Rack & pinion adjustment
Filter	Blue,Green,Yellow Filter
Light Source	Kholer Illumination 6V 20W halogen lamp with brightness control
Collector	For halogen lamp



A19.0209







**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China