



## 0.67x - 5x Trinocular Nikon SMZ745T Zoom Stereo Microscope with C-mount adapter A23.0710

Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A23.0710
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Magnification: 6.7x - 50x
- Applications: Industry|Gem|Medical|Clinic
- Certification: CE|Rohs
- Drawtube: Binocular|Trinocular
- Product Name: Stereo Zoom Microscope
- Theory: Stereo Microscope
- Highlight: **stereo dissecting microscope, stereo binocular microscopes**



### More Images



## Product Description

### 0.67x - 5x Trinocular Nikon SMZ745T Zoom Stereo Microscope with C-mount adapter A23.0710

#### Description:

- 0.67~5x Zoom Lens, with WF10x/22mm Eyepiece,
- Total Optical Magnification 6.7~50x
- Track Stand, No Light Source

A23.0710 Nikon SMZ745T Zoom Stereo Microscope Specification		T	TC
Head	Trinocular Head, 45° inclined, Interpupillary Adjustment 52-75mm		
Eyepiece	CW 10xB(FN 22)		
Zoom Lens	0.67 - 5X		
Zoom Ratio	1 : 7.5		
Magnification	1x-75x, Up to 3.35 - 300X With Optional Eyepieces and Auxiliary Objectives		
Working Distance	115mm		
C-mount	Built in a 0.55x C-mount Adapter		
Light Source	No Light Source		
Stand	Nikon Original Stand and Base	o	
	China-made Track Stand and Base		o



A23.0710





Back view of A23.0710



Splitting Ratio R:T=100:0 or 0:100

Zoom Lens: 0.67 - 5X



CW 10xB/(FN 22) Eyepiece

45° inclined, Interpupillary Adjustment  
52-75mm



**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China