



Nikon E200 Lab Wide Field Microscope CFI Optical System A12.0706

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A12.0706
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Laboratory
- Certification: CE
- Drawtube: Trinocular
- Illumination: LED
- Theory: Biological Microscope
- Highlight: binocular compound microscope,
optical compound microscope



More Images



Product Description

Nikon E200 Lab Wide Field Microscope CFI Optical System A12.0706

A12.0706 Nikon E200 Lab Microscope:


- Four magnification settings 40X, 100X, 400X, 1000X.
- CFI60 optics producing sharp and clear image.
- Ergonomic design brings comfortable operation
- Coaxial coarse/fine focusing, coarse motion torque adjustable, refocusing available

A12.0706 Biological Microscope		
Optical System	CFI Optical System	
Head	E2-TB Binocular Viewing Head, Inclined at 30°, Interpupillary 47-75mm	
Eyepiece	CFI W10×/20	
Objective	CFI E Plan Achromatic Objective 4x ,10x , 40x , 100x	
Nosepiece	Backward Quadruple Nosepiece	
Focusing	Coaxial Coarse and Fine Adjustment, Fine Division 0.002mm, Coarse Stroke 37.7mm per Rotation, Fine Stroke 0.2mm per Rotation	
Stage	Rectangle Stage 216×150mm, Moving Range 78×54mm	
Condenser	E2 Abbe Condenser NA1.25	
Illumination	6V/20W Halogen Lamp	
Optional Accessories		Item No.
Head	E2-TF Trinocular Viewing Head, Inclined at 30°, Interpupillary 47-75mm	
Phase Contrast	Phase Contrast Kit	
Attachment	Fluorescent Attachment	
	Dark Field Attachment	
	135 video Attachment	
	Photo Attachment	
	Digital video Attachment	




A12.0706 NIKON E200 Microscope





Optional parts:





Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China