



64 LED Light Microscope Accessory Blue / Red / Yellow / White A56.5204-04

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A56.5204-04
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Highlight: **compound microscope parts, parts of the microscope**

Product Description

64 LED Light Microscope Accessory Blue / Red / Yellow / White A56.5204-04

A56.5202-4 Color LED Ring Light, 60 LEDs

--64 LEDs Ring Light

--Four Colors Available : Blue, Red, Yellow, White, Four Color --Brightness Adjustable 9~100%

A56.5202-4 Four Color LED Ring Light	
Color	Blue, Red, Yellow, White, Four Color Separately Controlled By Switch
LED	64pcs
Inner Diameter	64mm
Light Intensity	60000Lum
Color Temperature	White Color 6400K
Power	3W
Input Voltage	Wide Range Voltage90-240V
Brightness Adjustment	0-100%
Focal Length	75mm
Working Distance	Effective Working Distance 40-160mm



A56.5202-4 Four Color LED Ring Light Source



Four Color Separately Controlled By Switch



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China