



Student Trinocular / Biogocial Compound Optical Microscope LED A11.1535

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1535
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Student
- Certification: CE
- Drawtube: Trinocular
- Illumination: LED
- Theory: Biological Microscope
- Manification: 40x~1000x
- Highlight: **binocular compound microscope, wide field microscope**



Product Description

Student Trinocular / Biogocial Compound Optical Microscope LED A11.1535

- 40x~1000x Magnification
- LED Light Source
- Coaxial Coarse & Fine Focusing Knob

A11.1535 Biological Microscope		A	B	C
Head	Monocular Head, Inclined at 45°, 360° Rotatable	o		
	Seidentopf Binocular Viewing Head, Inclined at 45°, Interpupillary 48-75mm		o	
	Seidentopf Trinocular Viewing Head, Inclined at 45°, Interpupillary 48-75mm			o
Eyepiece	WF10x/18mm			
Nosepiece	Quadruple Nosepiece			
Objective	Achromatic Objective 4X,10X,40X,100X			
Condensor	Abbe Condensor NA1.25 with Diaphragm, Filter			
Focusing	Coaxial Coarse and Fine Adjustment, Moving Range 22mm			
Working Stage	Double Layer Mechanical Rectangular Stage 140x130mm			
Light Source	LED Light, Brightness Adjustable			
Product Size	25*21.5*37cm			
Optional Accessories				
Eyepiece	WF16x/13mm			
	WF20x/11mm			
Nosepiece	Triple Nosepiece			
Light	6V20W Halogen Lamp, Brightness Adjustment			



A11.1535



Double Layers Mechanical Stage
140*130mm



A11.1535-A

A11.1535-B

A11.1535-C



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China