



A11.1002 Mechanical Stage LED Illumination Achromatic Binocular Biological Microscope

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1002
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month

OPTO-EDU



Product Specification

- Head: Binocular Head, Inclined 30°
- Eyepiece: Wide Field Eyepiece WF10x/18
- Objective: Achromatic Objective
- Nosepiece: Quadruple Nosepiece
- Focusing: Coaxial Coarse & Fine Focusing, Fine Division 0.004mm
- Stage: Double Layer Mechanical Stage
- Condenser: Abbe N.A.1.2 With Iris, Rack & Pinion Adjustment
- Illumination: S-LED Illumination, Brightness Adjustable
- Highlight: **Mechanical Stage Biological Microscope, Achromatic Binocular Biological Microscope, LED Illumination Biological Microscope**

OPTO-EDU



More Images

OPTO-EDU



OPTO-EDU



Product Description

Compact Design Student Biological Microscope
Monocular, Binocular Head Available, With WF10x/18mm Eyepiece
Quadruple Nosepiece With Achromatic Objective 4x10x40x100x
Attachable Mechanical Stage, Double Layer Mechanical Stage Optional
Coaxial Coarse & Fine Focusing, S-LED Illumination, Brightness Adjustable



OPTO-EDU (BEIJING) CO., LTD.

OPTO-EDU

F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China
Tel:+8610 88696020 Fax:+8610 88696085

A11.1002

Student Biological Microscope



A11.1002 Details

OPTO-EDU



Binocular Head, Inclined 30° for A11.1002-B

Monocular Head, Inclined 45°for A11.1002-A/C



Double Layer Mechanical Stage for A11.1002-B/C



Achromatic Objective 4x,10x,40x,100x



Abbe N.A.1.2 With Iris



S-LED Illumination, Brightness Adjustable

A11.1002 Specification

OPTO-EDU



A11.1002-A



A11.1002-B



A11.1002-C

A11.1002 Biological Microscope		A	B	C
Head	Monocular Head, Inclined 45°	•		•
	Binocular Head, Inclined 30°		•	
Eyeiece	Wide Field Eyeiece WF10x/18	•	•	•
Objective	Achromatic Objective 4x	•	•	•
	Achromatic Objective 10x	•	•	•
	Achromatic Objective 40x	•	•	•
	Achromatic Objective 100x	•	•	•
Nosepiece	Quadruple Nosepiece	•	•	•
Focusing	Coaxial Coarse & Fine Focusing, Fine Division 0.004mm, Coarse Stroke 37.7mm Per Rotation, Fine Stroke 0.4mm Per Rotation, Moving Range 24mm	•	•	•
Stage	Attachable Mechanical Stage, Size 110*120mm,Moving Range 60*30mm	•		
	Double Layer Mechanical Stage, Size 125*120mm,Moving Range 60*30mm		•	•

Condenser	Abbe NA1.2 With Iris, Spiral Up & Down	•		
	Abbe N.A.1.2 With Iris, Rack & Pinion Adjustment		•	•
Illumination	S-LED Illumination, Brightness Adjustable	•	•	•
Filter	Blue Filter	•	•	•
	Green Filter	•	•	•
Optional Accessories		Item No.		
Head	Duel Viewing Head	A53.1003-1		
Eyepiece	Wide Field Eyepiece P16x	A51.1002-16		
Illumination	Plan-concave Mirror	A56.1001-1		



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China