



LED Light Mono 1x Digital Optical Microscope USB output A32.5001

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A32.5001
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Laboratory
- Certification: CE
- Drawtube: Trinocular
- Illumination: LED
- Resolution: 2.0M
- Theory: Digital Optical Microscope
- Highlight: digital usb microscope,
optical microscope with digital camera



More Images



Product Description

LED Light Mono 1x Digital Optical Microscope USB output A32.5001

Description:

- 2.0M Camera Digital Head,
- USB2.0
- LED Light Source
- Widely apply to clinical diagnosis, teaching experiment, pathological test and other micro-fields.

A32.5001 Digital Mono Microscope, 1x

Resolution	2.0M
Main Body	1.0x
Auxiliary Objective	1.5x and 2.0x
Stage	185*120*28mm
Working Distance	120-55mm
Output	USB
Light Source	LED Light Source





Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China