



Monocular 40x - 400x Polarizing Light Microscope LED Light A15.0011

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: OPTO-EDU
- Certification: CE
- Model Number: A15.0011
- Minimum Order Quantity: 1
- Price: Negotiation
- Packaging Details: Carton
- Delivery Time: 1month
- Payment Terms: T/T, West Union and Paypal
- Supply Ability: 10000



Product Specification

- Application: Educational
- Certification: CE
- Drawtube: Monocular
- Illumination: LED Lamp
- Nosepiece: Triple
- Stage: Round Stage
- Highlight: student compound microscope, compound light microscope



Product Description

Monocular 40x - 400x Polarizing Light Microscope LED Light A15.0011

A15.0011 Polarizing Microscope, newly developed in 2015, with modern design and high optical quality, is ideal tools in educational field.

- Monocular Head
- Triple Nosepiece
- Round Double Layer Mechanical Working Stage
- N.A.1.25 Condenser, Iris Diaphragm
- Polarizer & Analyzer, Built-in Bertrand Lens, Compensator
- Colorful Focusing Knob, Light Collector and Power Switch
- LED Power On/Off Indicator
- LED Brightness Adjustment Indicator
- Cordless Rechargeable Battery Built-in

Specification

A15.0011 Polarizing Microscope	
Head	Monocular Head, 30° Inclined, 360° Rotatable
Eyepiece	WF10x/18mm, WF15x/13mm
Nosepiece	Triple Nosepiece, Ball Bearing Inner Locating
Objective	Achromatic Objective 4x, 10x, 40x (S)
Working Stage	Round, Diameter 142mm, Rotatable 360°, Division 1°, Vernier division 6, Lockable
Focusing	Coaxial Coarse & Fine Focus Knob, With Tension Adjustable And Limit Stopper, Minimum Division 0.002mm
Condenser	N.A. 1.25 Condenser, Iris Diaphragm, Filter
Polarizing	Analyzer : 0-90° Rotatable, Can Be Moved Out Of Optical Path For Single Polarizing Observe Bertrand Lens Built-in
	Compensator : λ Slip (First Class Red) , $1/4\lambda$ Slip, Quartz Wedge Plate (Optional)
Light Source	Polarizer : Mounted Under Abbe Condenser, 360° Rotatable
	High Light LED, Brightness Adjustable, With LED Indicator
Power	Cordless Built-in Rechargeable Battery
	100~240V Wide Range Voltage Power On/Off LED Indicator
Packing	Net Weight 4.80Kgs/pc,
	Packing 2 pcs/Carton, G.W. 10Kgs, Size 34*48*43cm



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China