



OPTO-EDU A59.5141 4K 8.0M HDMI USB Digital Camera



Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A59.5141
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month

OPTO-EDU



Product Specification

- Resolution: 8.0M 4K
- Sensor: 1/2.8" Sony CMOS, Pixel Size 1.45um
- System: Built-in System, Mouse Control To Take Photo, Record Video
- Setting: White Balance, Brightness, Exposure, Saturation, Contrast, Sharpness, Etc.
- Language: Chinese, English, Etc.
- Infrared: Cut-off 380-680nm
- Highlight: **4K HDMI USB digital microscope camera ,
8.0M USB microscope digital camera ,
HDMI microscope camera with 4K resolution**

for more products please visit us on cnoec.com

Product Description

8.0M 4K Digital Camera With 1/2.8" Sony CMOS Sensor
 All Kinds Output 4K HDMI+USB+TF Card +USB Disk
 HDMI Interface With Built-in Menu Controlled By Mouse For Setting & Measuring
 USB To PC Win11 Software StrangeView
 White Balance, Brightness, Exposure, Saturation, Contrast, Sharpness, etc.



A59.5141

4K HDMI+USB Digital Camera, 8.0M



A59.5141 Specification

OPTO-EDU



A59.5141

Features

- 8.0M 4K Digital Camera With 1/2.8" Sony CMOS Sensor
- All Kinds Output 2.4G/5G 4K HDMI+USB+TF Card +USB Disk
- HDMI Interface With Built-in Menu Controlled By Mouse For Setting & Measuring
- 2.4G/5G To IOS/Android App T-Cam, USB To PC Win11 Software StrangeView
- White Balance, Brightness, Exposure, Saturation, Contrast, Sharpness, etc.

Features

- 2.8" LCD Touch Screen 8.0M 4K Digital Camera With 1/1.8" Sony CMOS Sensor
- All Kinds Output 2.4G/5G WIFI+4K HDMI+USB+TF Card +USB Disk
- HDMI Interface With Built-in Menu Controlled By Mouse For Setting & Measuring
- 2.4G/5G WIFI To IOS/Android App T-Cam, USB To PC Win11 Software StrangeView
- White Balance, Brightness, Exposure, Saturation, Contrast, Sharpness, etc.



A59.5142

WIFI+HDMI+USB Digital Camera, 8.0M		A59.5141	A59.5142
Resolution	8.0M 4K	•	•
Sensor	1/2.8" Sony CMOS, Pixel Size 1.45um	•	
	1/1.8" Sony CMOS, Pixel Size 2.0um		•
LCD Screen	2.8" IPS Touch Screen, With Built-in System, Support Zoom In/Out, Take Picture, Record Video.	•	•
System	Built-in System, Mouse Control To Take Photo, Record Video, With All Setting, Measuring, Calibration Function	•	•
Setting	White Balance, Brightness, Exposure, Saturation, Contrast, Sharpness, etc.	•	•
Export	2.4G/5G WIFI To IOS/Android App T-Cam	•	•
	4K HDMI To Monitor	•	•
	Type C USB-UVC To Computer Measuring Software StrangeView	•	•
	USBx2 For Mouse & Keyboard Control, For USB Disk Storage	•	•
	TF Card Slot	•	•
Language	Chinese, English, etc.	•	•
Infrared	Cut-off 380-680nm	•	•
HDR	Wide Dynamic	•	•
Interface	C-Mount	•	•
Power	AC 95~265V, DC Adapter 5V2A	•	•



Opto-Edu (Beijing) Co., Ltd.

0086 13911110627

sale@optoedu.com

cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China