



OPTO-EDU A59.5142 8.0M 4K Digital Microscope Camera with 1/1.8" Sony CMOS Sensor and 2.8" LCD Touch Screen

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A59.5142
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Resolution: 8.0M 4K
- Sensor: 1/1.8" Sony CMOS, Pixel Size 2.0um
- System: Built-in System, Mouse Control To Take Photo, Record Video
- Setting: White Balance, Brightness, Exposure, Saturation, Contrast, Sharpness, Etc.
- Language: Chinese, English, Etc.
- Infrared: Cut-off 380-680nm
- Highlight: **8.0M 4K Resolution Digital Microscope Camera, 1/1.8 Sony CMOS Sensor Microscope Camera, 2.8 LCD Touch Screen HDMI Microscope Camera**

Product Description



Opto-Edu (Beijing) Co., Ltd.

2.8" LCD Touch Screen 8.0M
4K Digital Camera With 1/1.8"
Sony CMOS Sensor
All Kinds Output 2.4G/5G
WIFI+4K HDMI+USB+TF Card
+USB Disk
HDMI Interface With Built-in
Menu Controlled By Mouse For
Setting & Measuring
2.4G/5G WIFI To IOS/Android
App T-Cam, USB To PC Win11
Software StrangeView
White Balance, Brightness,
Exposure, Saturation, Contrast,
Sharpness, etc.



0086 13911110627

sale@optoedu.com

cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China

A59.5142

LCD Touch Screen 4K WIFI+HDMI+USB Digital Camera, 8.0M



A59.5142 Specification

OPTO-EDU



A59.5141

Features

- 8.0M 4K Digital Camera With 1/2.8" Sony CMOS Sensor
- All Kinds Output 2.4G/5G WIFI+4K HDMI+USB+TF Card +USB Disk
- HDMI Interface With Built-in Menu Controlled By Mouse For Setting & Measuring
- 2.4G/5G WIFI To IOS/Android App T-Cam, USB To PC Win11 Software StrangeView
- White Balance, Brightness, Exposure, Saturation, Contrast, Sharpness, etc.

Features

- 2.8" LCD Touch Screen 8.0M 4K Digital Camera With 1/1.8" Sony CMOS Sensor
- All Kinds Output 2.4G/5G WIFI+4K HDMI+USB+TF Card +USB Disk
- HDMI Interface With Built-in Menu Controlled By Mouse For Setting & Measuring
- 2.4G/5G WIFI To IOS/Android App T-Cam, USB To PC Win11 Software StrangeView
- White Balance, Brightness, Exposure, Saturation, Contrast, Sharpness, etc.



A59.5142

WIFI+HDMI+USB Digital Camera, 8.0M		A59.5141	A59.5142
Resolution	8.0M 4K	•	•
Sensor	1/2.8" Sony CMOS, Pixel Size 1.45um	•	
	1/1.8" Sony CMOS, Pixel Size 2.0um		•
LCD Screen	2.8" IPS Touch Screen, With Built-in System, Support Zoom In/Out, Take Picture, Record Video.	•	•
System	Built-in System, Mouse Control To Take Photo, Record Video, With All Setting, Measuring, Calibration Function	•	•
Setting	White Balance, Brightness, Exposure, Saturation, Contrast, Sharpness, etc.	•	•
Export	2.4G/5G WIFI To IOS/Android App T-Cam	•	•
	4K HDMI To Monitor	•	•
	Type C USB-UVC To Computer Measuring Software StrangeView	•	•
	USBx2 For Mouse & Keyboard Control, For USB Disk Storage	•	•
	TF Card Slot	•	•
Language	Chinese, English, etc.	•	•
Infrared	Cut-off 380-680nm	•	•
HDR	Wide Dynamic	•	•
Interface	C-Mount	•	•
Power	AC 95~265V, DC Adapter 5V2A	•	•



A59.5142 With A12.2263 Laboratory Biological Microscope

HDMI output sharp & clear image to large monitor, with built-in full function menu to set all parameters, calibrate, capture image, record video and measure, all controlled by mouse easily.

