



OPTO-EDU A59.5121 7"/9"/10.1" Button LCD Digital Camera 2.0M USB TF Card

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A59.5121
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Resolution: Photo 2M/5M/8M Video 1080p
- Sensor: 2.0M CMOS 1/2.8"
- Pixel Size: 2.9μm X 2.9μm
- Basic Function: Capture Photo, Record Video, Output To Computer + TF Card



Product Description

7"/9"/10.1" LCD Screen With Button Control, 2.0M CMOS 1/2.8"
 Type C TVC Output To Computer, With Measure Software StrangeView
 Photo Resolution 2M/5M/8M, Video Resolution 1080p, Capture & Saved In TF Card
 Including 18650 Rechargeable Battery (Not For Air Shipment) & 8G TF Card
 Optional Adapter: -E = Eyepiece Adapter (Dia.23.2mm), -C = C-Mount Adapter (1" 32UN Screw)



OPTO-EDU (BEIJING) CO., LTD.

OPTO-EDU

F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China
 Tel:+8610 88696020 Fax:+8610 88696085

A59.5121

7"/9"/10.1" Button LCD Digital Camera, 2.0M, USB+TF
 Card



A59.5121 Specification

OPTO-EDU




7"/9"/10.1" Button
 LCD Screen 2.0M

Output Type-C UVC
 + TF

With Built-in
 Measure Menu




Model	Button LCD Camera			2K Button LCD Camera			4K Touch LCD Camera	
	A59.5121-			A59.5122-			A59.5124-	
	7E/7C	9E/9C	10E/10C	7E/7C	9E/9C	10E/10C	7E/7C	10E/10C
LCD Size	7"	9"	10.1"	7"	9"	10.1"	7"	10.1"
LCD Resolution	1024*600						Touch Screen 1280*800	
Resolution	Photo 2M/5M/8M Video 1080p			Photo 2M/5M/8M/12M Video 1080p			Photo 5M/8M/12M Video 1080p	
Sensor	2.0M CMOS1/2.8"			2K 4.0M CMOS 1/3"			4K 8.0M CMOS 1/2.8"	
Pixel Size	2.9μm x 2.9μm			2μm x 2μm			1.45μm x 1.45μm	
Basic Function	Capture Photo, Record Video, Output To Computer +TF Card							
WIFI	-			2.4G, App T-Cam			2.4G/5G, App T-Cam	
HDMI	-			1080P-HDMI To Monitor			4K-HDMI To Monitor	
Mouse Control	-			Built-in Measure Menu			Built-in Adv. Measure Menu	
Type-C	Type-C Power In / USB-UVC Output To Computer, Measure Software Strange View							
USB	-			USBx1			USBx2	
TF Card	8G TF Card Support Max 64G						32G TF Card Support Max 64G	
Battery	18650-3.7V/2000MAH							
MINI USB	-			Connect AHD External Lens			-	
Adapter	-E = Eyepiece Adapter (Dia.23.2mm), -C = C-Mount Adapter (1" 32UN Screw)							
Power	Type-C Power In DC5V-2A							



-C = C-Mount Adapter
(1" 32UN Screw)



-E = Eyepiece Adapter
(Dia.23.2mm),



TF Card



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China