

OPTO-EDU A59.5122 Microscope Digital Camera with 2K Resolution 4.0M CMOS Sensor and 7/9/10.1 Inch LCD Screen

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A59.5122
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Resolution: Photo 2M/5M/8M/12M Video 1080p
- Sensor: 2K 4.0M CMOS 1/3"
- Pixel Size: 2μm X 2μm
- Basic Function: Capture Photo, Record Video, Output To Computer + TF Card
- Highlight: **2K Resolution Microscope Digital Camera, 4.0M CMOS Sensor LCD Digital Camera, 7/9/10.1 Inch LCD Screen Microscope Camera**



Product Description

OPTO-EDU A59.5122 7"/9"/10.1" 2K Button LCD Digital Camera

Key Features

- 7"/9"/10.1" 2K LCD Screen with Button Control & Mouse Control
- 4.0M CMOS 1/3" sensor for high-quality imaging
- Multiple output options: Type-C UVC, HDMI, 2.4G WIFI, and TF card
- PC Software StrangeView and Android/iOS APP T-Cam compatibility
- Photo resolutions: 2M/5M/8M; Video resolution: 1080p
- Includes 18650 rechargeable battery (not for air shipment) & 8G TF card
- Optional adapters: Eyepiece Adapter (-E) or C-Mount Adapter (-C)

**OPTO-EDU (BEIJING) CO., LTD.**

F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China
Tel:+8610 88696020 Fax:+8610 88696085

A59.5122

**7"/9"/10.1" 2K Button LCD Digital Camera, 4.0M,
USB+HDMI+WIFI+TF Card**



A59.5122 Specification

OPTO-EDU



7"/9"/10.1" 2K Button
LCD Screen 4.0M

With Built-in
Measure Menu

Output Type-C
UVC+HDMI+2.4G
WIFI+TF

PC Software
Strange View

Android / IOS APP
T-Cam



Model Comparison

Model	Button LCD Camera	2K Button LCD Camera	4K Touch LCD Camera
Model Number	A59.5121-	A59.5122-	A59.5124-
Available Sizes	7E/7C	9E/9C, 10E/10C	7E/7C, 10E/10C
LCD Size	7"	9", 10.1"	7", 10.1"
LCD Resolution	1024*600	Touch Screen 1280*800	Touch Screen 1280*800
Resolution	Photo: 2M/5M/8M Video: 1080p	Photo: 2M/5M/8M/12M Video: 1080p	Photo: 5M/8M/12M Video: 1080p
Sensor	2.0M CMOS1/2.8"	2K 4.0M CMOS 1/3"	4K 8.0M CMOS 1/2.8"
Pixel Size	2.9μm x 2.9μm	2μm x 2μm	1.45μm x 1.45μm
Basic Function	Capture Photo, Record Video, Output To Computer +TF Card		
WIFI	-	2.4G, App T-Cam	2.4G/5G, App T-Cam
HDMI	-	1080P-HDMI To Monitor	4K-HDMI To Monitor
Mouse Control	-	Built-in	Built-in
Measure Menu	Built-in	Adv. Measure Menu	Adv. Measure Menu
Type-C	Type-C Power In / USB-UVC Output To Computer, Measure Software Strange View		
USB	-	USBx1	USBx2
TF Card	8G TF Card Support Max 64G	32G TF Card Support Max 64G	32G TF Card Support Max 64G
Battery	18650-3.7V/2000MAH		
MINI USB	-	Connect AHD External Lens	-
Adapter	-E = Eyepiece Adapter (Dia.23.2mm), -C = C-Mount Adapter (1" 32UN Screw)		
Power	Type-C Power In DC5V-2A		



-C = C-Mount Adapter
(1" 32UN Screw)




-E = Eyepiece Adapter
(Dia. 23.2mm)



TF Card



Type-C To USB Computer
Win11 Software Strange View



4K HDMI With Built-in Menu
Mouse Control Setting & Measuring



iOS Android Google Play




2.4G/5G WIFI, App T Cam





A59.5122

HD 2K 9" Digital LCD Camera







Note: The included 18650 rechargeable battery is not approved for air shipment.



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China