



OPTO-EDU A59.5124 7"/10.1" 4K LCD Digital Camera 8.0M USB HDMI WIFI TF Card

Our Product Introduction

for more products please visit us on cnoec.com


Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A59.5124
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Resolution: Photo 5M/8M/12M Video 1080p
- Sensor: 4K 8.0M CMOS 1/2.8"
- Pixel Size: 1.45μm X 1.45μm
- Basic Function: Capture Photo, Record Video, Output To Computer +TF Card
- Highlight: **4K LCD digital microscope camera, USB HDMI microscope camera, WIFI TF card microscope camera**



Product Description

7"/10.1" 4K LCD Screen With Built-in Measure Menu, 8.0M CMOS 1/2.8"
 Output Type-C UVC+HDMI+2.4G WIFI+TF, PC Software StrangeView, Andorid/IOS APP T-Cam
 Photo Resolution 5M/8M/12M, Video Resolution 1080p, Capture & Saved In TF Card
 Including 18650 Rechargeable Battery (Not For Air Shipment) & 8G TF Card
 Optional Adapter: -E = Eyepiece Adapter (Dia.23.2mm), -C = C-Mount Adapter (1" 32UN Screw)



OPTO-EDU (BEIJING) CO., LTD.

OPTO-EDU

F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China
 Tel:+8610 88696020 Fax:+8610 88696085

A59.5124

7"/10.1" 4K LCD Digital Camera, 8.0M,
 USB+HDMI+WIFI+TF Card



A59.5124 Specification

OPTO-EDU



7"/10.1" 4K LCD
 Screen 8.0M

Output Type-C
 UVC+HDMI+2.4G
 WIFI+TF+USB3.0

With Built-in
 Measure Menu

PC Software
 Strange View

Android / IOS APP
 T-Cam



Button LCD Camera

2K Button LCD Camera

4K LCD Camera

Model	A59.5121-			A59.5122-			A59.5124-	
	7E/7C	9E/9C	10E/10C	7E/7C	9E/9C	10E/10C	7E/7C	10E/10C
LCD Size	7"	9"	10.1"	7"	9"	10.1"	7"	10.1"
LCD Resolution	1024*600						1280*800	
Resolution	Photo 2M/5M/8M Video 1080p			Photo 2M/5M/8M/12M Video 1080p			Photo 5M/8M/12M Video 1080p	
Sensor	2.0M CMOS/2.8"			2K 4.0M CMOS 1/3"			4K 8.0M CMOS 1/2.8"	
Pixel Size	2.9μm x 2.9μm			2μm x 2μm			1.45μm x 1.45μm	
Basic Function	Capture Photo, Record Video, Output To Computer +TF Card							
WIFI	-			2.4G, App T-Cam			2.4G/5G, App T-Cam	
HDMI	-			1080P-HDMI To Monitor			4K-HDMI To Monitor	
Mouse Control	-			Built-in Measure Menu			Built-in Adv. Measure Menu	
Type-C	Type-C Power In / USB-UVC Output To Computer, Measure Software Strange View							
USB	-			USBx1			USBx2	
TF Card	8G TF Card Support Max 64G						32G TF Card Support Max 64G	
Battery	18650-3.7V/2000MAH							
MINI USB	-			Connect AHD External Lens			-	
Adapter	-E = Eyepiece Adapter (Dia.23.2mm), -C = C-Mount Adapter (1" 32UN Screw)							
Power	Type-C Power In DC5V-2A							

A59.5124 Output

OPTO-EDU



-C = C-Mount Adapter
(1" 32UN Screw)




-E = Eyepiece Adapter
(Dia.23.2mm)



iOS

Android

Google Play



2.4G/5G WIFI, App T Cam




TF Card



USB 3.0 to U Disk



Type C TVC To Computer
Win11 Software Strange View



4K HDMI With Built-in Menu
Mouse Control Setting & Measuring



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China