



## OPTO-EDU A59.5124 7"/10.1" 4K Touch LCD Digital Camera 8.0M USB HDMI WIFI TF Card

### Our Product Introduction


#### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A59.5124
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



#### Product Specification

- Resolution: Photo 5M/8M/12M Video 1080p
- Sensor: 4K 8.0M CMOS 1/2.8"
- Pixel Size: 1.45μm X 1.45μm
- Basic Function: Capture Photo, Record Video, Output To Computer +TF Card



## Product Description

7"/10.1" 4K LCD Touch Screen With Built-in Measure Menu, 8.0M CMOS 1/2.8"  
 Output Type-C UVC+HDMI+2.4G WIFI+TF, PC Software StrangeView, Andorid/IOS APP T-Cam  
 Photo Resolution 5M/8M/12M, Video Resolution 1080p, Capture & Saved In TF Card  
 Including 18650 Rechargeable Battery (Not For Air Shipment) & 8G TF Card  
 Optional Adapter: -E = Eyepiece Adapter (Dia.23.2mm), -C = C-Mount Adapter (1" 32UN Screw)



OPTO-EDU (BEIJING) CO., LTD.

OPTO-EDU

F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China  
 Tel:+8610 88696020 Fax:+8610 88696085

## A59.5124

7"/10.1" 4K Touch LCD Digital Camera, 8.0M,  
 USB+HDMI+WIFI+TF Card



## A59.5124 Specification

OPTO-EDU



7"/10.1" 4K LCD Touch  
 Screen 8.0M

Output Type-C  
 UVC+HDMI+2.4G  
 WIFI+TF+USB3.0

With Built-in  
 Measure Menu

PC Software  
 Strange View

Android / IOS APP  
 T-Cam



|       | Button LCD Camera | 2K Button LCD Camera | 4K Touch LCD Camera |
|-------|-------------------|----------------------|---------------------|
| Model | A59.5121-         | A59.5122-            | A59.5124-           |

[illegible]

## A59.5124 Output

OPTO-EDU



-C = C-Mount Adapter  
(1" 32UN Screw)




-E = Eyepiece Adapter  
(Dia.23.2mm)



iOS

Android

Google Play



2.4G/5G WIFI, App T Cam




TF Card



USB 3.0 to U Disk



Type C TVC To Computer  
Win11 Software Strange View



4K HDMI With Built-in Menu  
Mouse Control Setting & Measuring



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China