



OPTO-EDU A22.1508 2x 4x Binocular LED Pole Stand Step Zoom Stereo Microscope

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A22.1508
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month

OPTO-EDU



Product Specification

- Head: Binocular Head, 45° Inclined, Interpupillary Distance 55~75mm
- Eyepiece: WF10x, 1 Pair, A51.1523-10
- Objective: 2x/4x
- Focus: Coarse Focusing Knob, Focusing Range 40mm
- W.D.: Working Distance 57mm
- Stage: Frosted Glass Plate, Black/White Plate, Dia.60mm
- Light Source: Upper/Bottom LED, Brightness Adjustable Separately
- Power: AC 110-220V

OPTO-EDU



More Images

OPTO-EDU



OPTO-EDU



OPTO-EDU



Product Description

Step Zoom Stereo Microscope 1x3x, 1x2x, 2x4x Optional
Binocular Head Inclined 45° Interpupillary Distance Adjustable
WF10x Eyepiece 1 Pair, With Objective Reach Magnification 10x-40x
Compact Designed Pole Stand With Coarse Focusing Knob
Transmit & Reflect LED Light Source Brightness Adjustable



OPTO-EDU (BEIJING) CO., LTD.

OPTO-EDU

F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China
Tel: +8610 88696020 Fax: +8610 88696085

A22.1508

Stereo Microscope, 2x/4x



A22.1508 Specification

OPTO-EDU




WF10x Eyepiece 1 Pair, With Objective
Reach Magnification 10x-40x

Binocular Head Inclined 45° Interpupillary
Distance Adjustable

Transmit & Reflect LED Light Source
Brightness Adjustable

Compact Designed Pole Stand With
Coarse Focusing Knob

A22.1508 Stereo Microscope		
Head	Binocular Head, 45° Inclined, Interpupillary Distance 55~75mm, Left Eyepiece Tube Diopter Adjustable +/-5°	●
Eyepiece	WF10x, 1 Pair, A51.1523-10	●
	WF20x, 1 Pair, A51.1523-20	○
Objective	2x/4x	●
	1x/2x	○
	1x/3x	○
Focus	Coarse Focusing Knob, Focusing Range 40mm	●
W.D.	Working Distance 57mm	●
Stage	Frosted Glass Plate, Black/White Plate, Dia.60mm	●
Light Source	Upper/Bottom LED, Brightness Adjustable Separately	●
Power	AC 110-220V	●
Note: ● Means Standard Outfit, ○ Means Optional		





Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China