



OPTO-EDU A36.3602 Professional 110X LCD Digital Measuring Microscope 11.6 inch 5MP CMOS Camera Stereo Video Microscope

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A36.3602
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month

OPTO-EDU



Product Specification

- Highlight: **Professional LCD digital measuring microscope , 11.6 inch stereo video microscope , 5MP CMOS camera microscope**

for more products please visit us on cnoec.com

Product Description

OPTO-EDU A36.3602 Professional 110X LCD Digital Measuring Microscope

11.6 inch 5MP CMOS Camera Stereo Video Microscope

The A36.3602 LCD stereo microscope integrates a zoom optical module with high-definition camera, designed for convenient use in narrow and small workspaces. Featuring ergonomic internal optical corner design and adjustable display screen angle, it allows operators to view the display comfortably. Through the HDMI HD interface, it can connect to an additional display simultaneously.

A36.3602 LCD Digital Stereo Microscope Specifications

Magnification	17X-110X
Objectives	0.75X-4.5X stereo zoom, zoom ratio 6.4:1
Display	11.6 Inch, 1920H*1080V, IPS screen, detachable
Imaging sensor	1/2.8 Inch SONY CMOS
Real-time image frame rate	60fps@1920*1080
Photograph	5 Mega-pixel
Recording	45fps@1920*1080
Output Interface	HDMI
Input Interface	Double USB2.0
Camera Parameters	Exposure, White-balance, Gain, Gamma, Brightness, Contrast, Hue, Saturation, Gray
OSD Operation	Mouse
Storage Function	Snap, Record, USB Storage
Measurement	Point, Line, Parallel lines, Broken lines, Angle, Diameter circle, Radius circle, Three-point circle, Distance from circle to line, Concentric circles, Double circle, Vertical line, Arbitrary line, Polygon, Arc, Mark, Screenshot, Table data export, Calibration, Parameter Settings
Other Functions	Flip horizontally, Flip vertically, Freeze
Scale Ruler	3 kinds, customizable colors, line widths, display switches
Grid Lines	8 group, customizable colors, line widths, display switches
Languages	Simplified Chinese, Traditional Chinese, English, French, Russian, Japanese, Italian, German
Auxiliary Lens	0.5x, 2.0x
Light Source	56-LED ring light, brightness adjustable
Base	Various bases available according to customer requirements

A36.3602 LCD Stereo Microscope With Menu To Take Photo, Video. Adjust Parameter, and Measure, Controlled By Mouse.



Microscope Optical Data Sheet

Objective	Standard	Auxiliary 0.5X	Auxiliary 2X
W.D.	100mm	165mm	30mm
11.6" Magnification	17.2x-110.5x	8.6x-55.2x	34.4x-221x
View Field	18.5x10.5-2.8x1.5	37.5x21-5.5x3.0	9.5x5.5-1.5x0.7

The coin under the stereo microscope:



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China