



## OPTO-EDU A59.4981 15.6" 4K LCD Monitor For Microscope

Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A33.1550
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month

OPTO-EDU



### Product Specification

- Highlight: 15.6 inch 4K LCD microscope monitor , digital LCD microscope with 4K resolution , OPTO-EDU A59.4981 microscope display

OPTO-EDU



### More Images

OPTO-EDU



OPTO-EDU






## Grid Measurement

Draw customizable grids with adjustable color and line width to assist measurement work.

## Auto Focus

Built-in motor supports real-time Auto Focus with three modes: Continuous, One Time, and Manual.



## Illumination Controller


Software-controlled illumination angle and intensity for optimal viewing conditions.

## EDF Function

Extended Depth of Field function stacks multiple image layers to create fully focused, clear images.

## Large Image Stitching

Supports Auto Stitching, Manual Stitching, and Extend Column Stitching to create large composite images.



### Combined Color Images

Merge different fluorescent bands images into a single composite image for comprehensive analysis.



**Opto-Edu (Beijing) Co., Ltd.**



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China