



OPTO-EDU A59.3506 4K HDMI+USB Digital Camera, 8M, Mouse Measure+USB Measure

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A59.3506
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month

OPTO-EDU



Product Specification

- Highlight: 4K HDMI digital microscope camera, 8MP USB microscope camera, mouse measure microscope camera

OPTO-EDU



More Images

OPTO-EDU



OPTO-EDU



Product Description

OPTO-EDU A59.3506 4K HDMI+USB Digital Camera, 8M, Mouse Measure+USB Measure

Product Overview

- 4K 8MP resolution with 1/1.8" Sony IMX678 CMOS sensor
- Triple output capability: HDMI + USB + USB Disk
- HDMI interface with built-in user interface menu supporting USB mouse control
- HDMI software supports photo & video capture, measurement, and calibration
- Ratio setting: Input objective, eyepiece, and monitor size to show real-time magnification
- USB UVC output to computer with S-Eye Measure software for Windows



A59.3506

4K 8M HDMI+USB Digital Camera






Technical Specifications

Items	A59.3506 4K 8M HDMI+USB Digital Camera
Camera	8.0M 4K Digital Camera
Sensor	1/1.8" Sony IMX678 CMOS
Resolution	3840x2160@30fps
Pixel Size	2umx2um
Exposure Capability	Auto or Manual
Exposure Time	1~999ms
White Balance	Auto, One-Push, or Manual RGB Adjustment
Power Supply	Power Adaptor DC 12V 1A
Digital Zoom	1~10x
Data Output	HDMI Output x1 to Monitor 4K @ 30fps USB 2.0 Output x3 to: -- PC 4K @ 30fps -- USB Disk -- Connect USB Mouse Support USB+HDMI Sync Output
Software	[HDMI Output] Built-in Software Support Mouse Control, With NO NEED For Computer: -- Ratio Setting (Automatic Calculation of Magnification) -- Photo, Video Capture & Store (.jpg or .avi) -- All Kinds Measure, Mark, Analysis Functions -- Mark & Measure Results Capture & Store (.jpg or .xlsx) -- Calibration -- Auto Edge Detection -- Settings [USB Output] S-EYE 2.0 Support Windows System -- Support All Above Functions Language Supported: English, Simplified Chinese, Traditional Chinese, Japanese, Korean
Interface	Standard C-Mount Port
Packing List	-- Digital Camera x1 -- 12V1A Power Adaptor x1 -- HDMI Cable x1 -- USB Cable x1 -- USB Mouse x1

Product Size & Weight	6(W)x6(L)x6.7(H)cm, 0.2Kg
Package Size & Weight	17.5x19x9cm, 0.5Kg

A59.3506 Structure

OPTO-EDU





♦ **HDMI output to TV monitor**


Easy USB Mouse Control, USB flash disk storage, Max frame rate: 30fps

♦ **USB output to Computer**


Camera support USB output, it can work with computer, Max frame rate: 30fps We will provide S-EYE2.0 software for buyers (S-EYE2.0 only support windows system)




◆ USB UVC Output To Computer With S-Eye Measure Software for Windows



◆ HDMI Software Support Photo & Video Capture, Measure, Calibration





Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China