



7x - 45x Binocular Stereo Optical Microscope with Track Stand LED Light A23.0906 - J3L

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A23.0906 - J3L
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Magnification: 7x-45x
- Applications: Industry|Gem|Medical|Clinic
- Certification: CE|Rohs
- Stand: Track Stand
- Product Name: Stereo Zoom Microscope
- Theory: Stereo Microscope
- Highlight: **stereo zoom microscopes,
stereo dissecting microscope**

Product Description

7x - 45x Binocular Stereo Optical Microscope with Track Stand LED Light A23.0906 - J3L

- Ergonomic design can make users feel comfortable after long-time operation
- Standard zoom magnification 7X-45 with zoom ratio 1:6.5. Standard working distance is 100mm
- Configured with WF10X/20mm eyepiece.
- Full set of eyepiece, auxiliary objectives and other accessories can fulfill different needs.

A23.0906-J3L Zoom Stereo Microscope, 0.7x-4.5x Specification		B	T
Head	Binocular, 45° Inclined, 360° Rotatable, +/-5° Diopter Adjustable, Interpupillary Adjustment 54~75mm Trinocular, 45° Inclined, 360° Rotatable, +/-5° Diopter Adjustable, Interpupillary Adjustment 54~75mm	o	o
Zoom Range	0.7~4.5x		
Zoom Ratio	1:6.5		
Eyepiece	WF10x/20mm, Both with Eye-cup		
Magnification	7~45x, Up to 3.5~180x with Auxiliary Lens		
Working Distance	100mm, With Standard 1x Lens		
Field of View	1~95.2mm		
Stand	Big Fan Track Stand, With Focusing Knob		
Illumination	Up & Bottom High Brightness 3W LED Light Source		
Optional Accessories			
Eyepiece	WF15x/15mm	A51.0922-15	
	WF20x/10mm	A51.0922-20	
	0.5x	A52.0927-05	
	0.75x	A52.0927-75	
Auxiliary Lens	1.5x	A52.0927-15	
	2.0x	A52.0927-20	



A23.0906-J3L-B

www.cnoec.com



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China