



Binocular Stereo Optical Microscope With Pole Stand A23.0903 - D

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A23.0903 - D
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Magnification: 6.7x-45x
- Applications: Industry|Gem|Medical|Clinic
- Certification: CE|Rohs
- Drawtube: Binocular|Trinocular
- Product Name: Stereo Zoom Microscope
- Theory: Stereo Microscope
- Highlight: **stereo dissecting microscope, stereo binocular microscopes**



More Images



Product Description

Binocular Stereo Optical Microscope With Pole Stand A23.0903 - D

- Ergonomic design can make users feel comfortable after long-time operation
- Standard zoom magnification 6.7 X-45 with zoom ratio 1:6.5. Standard working distance is 100mm
- Configured with WF10X/22mm eyepiece.
- Full set of eyepiece, auxiliary objectives and other accessories can fulfill different needs.

A23.0903-D Zoom Stereo Microscope, 0.67x-4.5x Specification		BTB1 T1
Head	Binocular, 45° Inclined, 360° Rotatable, +/-5° Diopter Adjustable, Interpupillary Adjustment 54~75mm	o o
	Trinocular, 45° Inclined, 360° Rotatable, +/-5° Diopter Adjustable, Interpupillary Adjustment 54~75mm	o o
Zoom Range	0.67~4.5x	
Zoom Ratio	1:6.5	
Eyepiece	WF10x/20mm, Both with Eye-cup	
Magnification	6.7~45x,	
Working Distance	105mm	
Stand	Pole Stand, With Focusing Knob, Focusing Range 50mm	o o
	Track Stand, With Focusing Knob, Focusing Range 50mm	o o
Illumination	No Light Source	
A23.0903-D- Zoom Stereo Microscope, 0.67x-4.5x Optional Accessories		
Eyepiece	WF15x/15mm	A51.0923-15
	WF20x/10mm (Micrometer is available, Optional)	A51.0923-20
Auxiliary Lens	0.5x, W.D. 165mm	A52.0928-05
	2x, W.D. 30mm	A52.0928-75
Adapter	0.5x, C-Mount	A55.3601-05
	1x, C-Mount	A55.3601-1



A23.0903-D-B Binocular Stereo Microscope



A23.0903-D-T Trinocular Stereo Microscope



A23.0903-D-B1 Binocular Stereo Microscope



A23.0903-D-T1 Trinocular Stereo Microscope



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China