



## 200 X USB Digital Optical Microscope for Andorid Mobile Phone A34.5011

### Our Product Introduction

#### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A34.5011
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



#### Product Specification

- Applications: Laboratory
- Certification: CE
- Magnification: 100x -200x
- Illumination: LED
- Resolution: 1600\*1200
- Theory: Digital Optical Microscope
- Highlight: **digital usb microscope,**  
**portable digital microscope**



## Product Description

### USB Polarizing Digital Optical Microscope 500x LED Light A34.5012

- Magnification:100-200x
- Operating System:Windows &, Mac (USB port) Andriod

| A34.5011 USB Digital Microscope For Andorid Mobile Phone |                                                                                                         |
|----------------------------------------------------------|---------------------------------------------------------------------------------------------------------|
| Magnification                                            | 10x~200x, (at 5mm on 17" screen by 640 x 480 resolution)                                                |
| Sensor                                                   | High-quality CMOS sensor, 2.0M pixels                                                                   |
| Resolution                                               | 1600 X 1200                                                                                             |
| Focal distance:                                          | 5mm - $\infty$                                                                                          |
| Interface                                                | USB2.0                                                                                                  |
| Frame rate:                                              | 30 frames/sec (CIF and VGA)                                                                             |
| Support System                                           | Windows &, Mac (USB port)<br>Andriod ( A standard accessory: Special OTG line and EHECAM.apk software ) |
| USB cable length:                                        | 1.5 meters                                                                                              |





**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ [sale@optoedu.com](mailto:sale@optoedu.com)

🌐 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China