



## Monocular Head Cordless Biological Microscope For Student A11.0021

### Our Product Introduction

#### Basic Information

- Place of Origin: Korea
- Brand Name: OSUN
- Certification: CE, Rohs
- Model Number: A11.0021
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



#### Product Specification

- Application: Student
- Head: Monocular
- Nosepiece: Triple
- Illumination: LED Lamp
- Power: Cordless
- Condensor: N.A. 1.25
- Highlight: **high powered microscope ,  
biological compound microscope**



## Product Description

### Monocular Head Cordless Biological Microscope For Student A11.0021

#### Made in Korea!

**A11.0021 Cordless Student Microscope**, newly developed in 2015, with modern design and high optical quality, is ideal tools in educational field.

- Monocular Head Rotatable
- Triple Nosepiece
- Double Layer Mechanical Working Stage
- N.A.1.25 Condenser, Iris Diaphragm
- Colorful Focusing Knob, Light Collector and Power Switch
- LED Power On/Off Indicator
- LED Brightness Adjustment Indicator
- Cordless Rechargeable Battery Built-in

A11.0021 Cordless Student Microscope		A11.0021	-60	-100
Head	Monocular Head, 30° Inclined, 360° Rotatable			
Eyeiece	WF10x/18mm Wide Field Eyeiece			
Nosepiece	Triple Nosepiece, Ball Bearing Inner Locating			
Objective	Achromatic Objective 4x, 10x, 40x (S)	o		
	Achromatic Objective 10x, 40x (S), 60x (S)		o	
	Achromatic Objective 10x, 40x (S), 100x (S)			o
Working Stage	Double Layer Mechanical Stage, Size 125*120mm, Moving Range 72*30mm			
Focusing	Coaxial Coarse & Fine Focus Knob, With Tension Adjustable And Limit Stopper, Minimum Division 0.002mm			
Condenser	N.A. 1.25 Condenser, Iris Diaphragm, Filter			
Light Source	High Light LED, Brightness Adjustable, With LED Indicator			
Power	Cordless Built-in Rechargeable Battery 100~240V Wide Range Voltage Power On/Off LED Indicator			
Packing	Net Weight 4.28 Kgs/pc, Packing 2 pcs/Carton, G.W. 10Kgs, Size 34*48*43cm			



A11.0021 Student Microscope Front View



A11.0021 Student Microscope Back View



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



[sale@optoedu.com](mailto:sale@optoedu.com)



[cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China