



## Portable Bonocular Polarizing Light Microscope Cordless 40x - 400x A15.0022

### Our Product Introduction

#### Basic Information

- Place of Origin: China
- Brand Name: OPTO-EDU
- Certification: CE
- Model Number: A15.0022
- Minimum Order Quantity: 1
- Price: Negotiation
- Packaging Details: Carton
- Delivery Time: 1 month
- Payment Terms: T/T, West Union and Paypal
- Supply Ability: 10000



#### Product Specification

- Application: Student
- Certification: CE
- Drawtube: Binocular
- Illumination: LED Lamp
- Magnification: 40x-100x
- Stage: Round Stage
- Highlight: student compound microscope, highest magnification microscope



## Product Description

### Portable Bonocular Polarizing Light Microscope Cordless 40x - 400x A15.0022

A15.0022 Polarizing Microscope, newly developed in 2015, with modern design and high optical quality, is ideal tools in educational field.

- Binocular Head
- Triple Nosepiece
- Round Double Layer Mechanical Working Stage
- N.A.1.25 Condenser, Iris Diaphragm
- Polarizer & Analyzer, Built-in Bertrand Lens, Compensator
- Colorful Focusing Knob, Light Collector and Power Switch
- LED Power On/Off Indicator
- LED Brightness Adjustment Indicator
- Cordless Rechargeable Battery Built-in

### Specification

<b>A15.0022 Polarizing Microscope</b>	
Head	Binocular Head, 30° Inclined, 360° Rotatable
Eyeiece	WF10x/18mm, WF15x/13mm
Nosepiece	Triple Nosepiece, Ball Bearing Inner Locating
Objective	Achromatic Objective 4x, 10x, 40x (S)
Working Stage	Round, Diameter 142mm, Rotatable 360°, Division 1°, Vernier division 6, Lockable
Focusing	Coaxial Coarse & Fine Focus Knob, With Tension Adjustable And Limit Stopper, Minimum Division 0.002mm
Condenser	N.A. 1.25 Condenser, Iris Diaphragm, Filter
Polarizing	<b>Analyzer :</b> 0-90° Rotatable, Can Be Moved Out Of Optical Path For Single Polarizing Observe
	Bertrand Lens Built-in
	Compensator : $\lambda$ Slip (First Class Red) , 1/4 $\lambda$ Slip, Quartz Wedge Plate
Light Source	<b>Polarizer :</b> Mounted Under Abbe Condenser, 360° Rotatable, With Scale, Lockable
	High Light LED, Brightness Adjustable, With LED Indicator
Power	Cordless Built-in Rechargeable Battery
	100~240V Wide Range Voltage
Packing	Power On/Off LED Indicator
	Net Weight 4.80Kgs/pc,
	Packing 2 pcs/Carton, G.W. 10Kgs, Size 34*48*43cm



**Opto-Edu (Beijing) Co., Ltd.**



0086 13911110627



[sale@optoedu.com](mailto:sale@optoedu.com)



[cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China