



5.0 M WIFI Digital Microscope Monocular 365X Magnification A34.4903 - B

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A34.4903 - B
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: WIFI
- Certification: CE|Rohs
- Drawtube: Monocular
- Illumination: LED
- Max Magnification: 365X
- Theory: Digital Optical Microscope
- Highlight: portable digital microscope,
optical microscope with digital camera



Product Description

5.0 M WIFI Digital Microscope Monocular 365X Magnification A34.4903 - B

Description:

A34.4903-B Digital Biological Microscopes, newly designed in 2016, is equipped with 5.0M CMOS digital head, achromatic objectives and wide field eyepieces. It can provide clear image, wide field, convenience operation, can be used for medicine agriculture area widely. They are the ideal instruments in medical treatment, teaching demonstration.

A34.4903 WIFI Miroscope	
Microscope System	--0.7~5x, Continus Zoom Lens
	--Objective Lens 1.0x
	--Working Distance 105mm
	-- Pole Stand, With Focus Arm
	--Standard C-Mount
Digital Camera	--No Light Source
	5.0M CMOS High Resolution Digital Camera
Sensor	1/2.5" CMOS, True Color
Resolution	2592*1944 Max.
Pixel Size	2.2x2.2μm
Frame Rate	@640X 480: 30fps; @2592X1944: 15fps 720P HD Video Up To 30fps
White Balance	Auto
Exposure	Manual/Auto
Image Format	JPG/BMP/PNG/GIF
Video Format	AVI
Output	USB2.0 / Wifi Wireless
Software	Capturing Software, Free Driver
System	IOS/Android/Microsoft Windows XP/2003/Vista/Win7/Win8
Lens Mount	Standard C-mount



A34.4903-B WIFI Miroscope





Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China