



Phase Contrast Microscope

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE,
- Model Number: A19.0209-40
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Certification: CE
- Drawtube: Binocular
- Illumination: Halogen Lamp
- Theory: Phase Contrast Microscope
- Condensor: N.A. 1.25
- Nosepiece: Backward Quadruple
- Highlight: **phase contrast light microscope, inverted phase contrast microscope**

Product Description

Phase Contrast Microscope with Halogen Lamp Illumination A19.0209-40

- Compact Design For Appearance And Structure
- Ergonomic Design, Easy Operation, Outstanding Optical Performance, Clear and Bright Image
- Specially Developed for Phase Contrast Observation, High Contrast, High Resolution, Highlight the Cell Outline and Internal Structure
- High Precision Focusing Device and Convenient Illumination Adjustment

A19.0209-40 Phase Contrast Microscope	
Magnification	40x~1000x
Head	Sliding Binocular Head, 30° Inclined
Eyepiece	Wide Field WF10x/18mm
Nosepiece	Backward Quadruple Nosepiece
Objective	Achromatic 4x/0.10, 10x/0.25, 100xs/1.25(Oil)
Stage	Double Layer Mechanical Stage, 140*140mm, Mvoing Range 50*75mm
Focusing	Coaxial Coarse & Fine Focusing. Fine Focusing Scale 0.002mm, Coarse Focusing Stroke 37.7mm/Rotation, With Tension Adjustment, With Up Limit, Base Bearing Black Plastic Cap
Condenser	N.A.1.25 Abbe Condenser, Iris Diaphargm, With Filter Holder
Filter	Blue, Yellow, Green
Phase Contrast	Phase Contrast Objective 40x Phase Contrast Slide 40x Centering Telescope 10x
Light Source	Built-in Illumination 6V20W Halogen Lamp, 100~240V



A19.0209-40 Phase Contrast Microscope



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China