



Student Biological Microscope With Blue & Green Filter A11.1002

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1002
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Product Name: Microscope
- Theory: Biological Microscope
- Focusing: Coaxial Coarse & Fine Focusing
- Applications: Student; Laboratory
- Nosepiece: Quadruple Nosepiece
- Magnification: 40x-1000x
- Highlight: **high powered microscope ,
biological compound microscope**



More Images



Product Description

Specifications

A11.1002 Biological Microscope		A	B	C
Head	Monocular Head, Inclined 45°	●		●
	Binocular Head, Inclined 30°		●	
Eyepiece	Wide Field Eyepiece WF10x/18	●	●	●
Objective	Achromatic Objective 4x	●	●	●
	Achromatic Objective 10x	●	●	●
	Achromatic Objective 40x	●	●	●
	Achromatic Objective 100x	●	●	●
Nosepiece	Quadruple Nosepiece	●	●	●
Focusing	Coaxial Coarse & Fine Focusing, Fine Division 0.004mm, Coarse Stroke 37.7mm Per Rotation, Fine Stroke 0.4mm Per Rotation, Moving Range 24mm	●	●	●
Stage	Attachable Mechanical Stage, Size 110*120mm, Moving Range 60*30mm	●		
	Double Layer Mechanical Stage, Size 125*120mm, Moving Range 60*30mm		●	●
Condenser	Abbe NA1.2 With Iris, Spiral Up & Down	●		
	Abbe N.A.1.2 With Iris, Rack & Pinion Adjustment		●	●
Illumination	S-LED Illumination, Brightness Adjustable	●	●	●
Filter	Blue Filter	●	●	●
	Green Filter	●	●	●
Optional Accessories		Item No.		
Head	Dual Viewing Head	A53.1003-1		
Eyepiece	Wide Field Eyepiece P16x	A51.1002-16		
Illumination	Plan-concave Mirror	A56.1001-1		

A11.1002 Biological Microscope



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

