




## OPTO-EDU A32.6401 Digital Stereo Microscope 0.7x~4.5x 12M HDMI

### Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)


#### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A32.6401
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month





#### Product Specification

- Camera Pixel: 12MP
- Video Resolution: 1080p/60fps
- Data Output: HDMI
- WDR: Yes
- Storage: SD Card Storage (Max 32G)
- External Interface: Projector/Monitor With HDMI Port
- Highlight: **usb digital microscope,**  
**infinity plan microscope**



#### More Images



## Product Description

0.7x~4.5x Zoom Lens 14x-90x Total Magnification On 10.6" LCD

Full HD 12.0M Image And 1080p/30fps Video

Upper + Side LED Light Brightness Adjustable

SD Card Slot And HDMI Output To Monitor Available

Black Pole Stand 38\*25cm

A32.6401, A36.6401 Series Digital Microscope			A32.6401 -12M	A36.6401 -12M	A36.6401 -12MH
LCD Screen	Size	10.6" IPS LCD		●	●
	Resolution	1366*768		●	●
Digital Camera	Camera Pixel	12MP	●	●	●
	Video Resolution	1080p/60fps	●	●	●
	Data Output	HDMI	●		●
	WDR	Yes	●	●	●
	Function	Auto Exposure, Exposure Offset Auto White Balance Photo, Video Recording	●	●	●
	Storage	SD Card Storage (Max 32G)	●	●	●
	External Interface	Projector/Monitor With HDMI Port	●		●
Optical	Objective	0.7x~4.5x Zoom Lens	●	●	●
	Total Magnification	14x~90x On 10.6" LCD Screen	●	●	●
	Illumination	Side/Up LED Light, Brightness Adjustable	●	●	●
	Focus Range	60mm	●	●	●
	Working Distance	88mm	●	●	●
Stand	Column	Round Pole, Dia.25mm	●	●	●
	Stand	380mm*250mm*20mm, Black 320mm*260mm*20mm, White (Same Cost)	●	●	○
Optional	A52.6401-0.5	Auxiliary Lens 0.5x	○	○	○
	A52.6401-2.0	Auxiliary Lens 2.0x	○	○	○
	Boot Logo		○	○	○
	WiFi SD Card, Upgrade To WiFi Microscope		○	○	○
	SD Card		○	○	○
Power	DC 12V1A, Wide Range 100~240V AC Adapter		●	●	●

Remark: "●" Means Standard Outfit, "○" Means Optional



**0.7x~4.5x Zoom Lens 14x-90x  
Total Magnification On 10.6" LCD**

**Full HD 12.0M Image  
And 1080p/30fps Video**



## Product Details



## Product Size And Application





## Product Specification



A32.6401-12M



A36.6401-12M



A32.6401-12MH



**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ [sale@optoedu.com](mailto:sale@optoedu.com)

🌐 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China