



Monocular Student Biological Microscope With Separate Coarse & Fine Focusing A11.1003

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1003
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Product Name: Microscope
- Theory: Biological Microscope
- Focusing: Separate Coarse & Fine Focusing
- Applications: Student; Laboratory
- Head: Monocular Head
- Magnification: 40x-1000x
- Highlight: **high powered microscope ,
biological compound microscope**



More Images



Product Description

Specifications

A11.1003 Biological Microscope		A	B	C
Head	Monocular Head, Inclined 45°	•	•	•
Eyeiece	Wide Field Eyepiece WF10x/18	•	•	•
Objective	Achromatic Objective 4x	•	•	•
	Achromatic Objective 10x	•	•	•
	Achromatic Objective 40x	•	•	•
	Achromatic Objective 100x			•
Nosepiece	Triple Nosepiece	•	•	
	Quadruple Nosepiece			•
Stage	Plain Stage With Slide Clips 120*130mm	•	•	
	Attachable Mechanical Stage 110*120mm/ 60*30mm			•
Focusing	Separate Coarse & Fine Focusing	•	•	•
Condenser	Single Lens With Disc Diaphragm	•	•	
	Abbe NA 1.2 With Iris Diaphragm And Filter			•
Illumination	Plan-concave Mirror	•		
	Built-in Illumination, Incandescent Lamp 110V/ 220V/ 15W		•	•
Optional Accessories		Item No.		
Eyeiece	Wide Field Eyepiece P16x	A51.1002-16		

A11.1003 Biological Microscope



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China