



Inner And Out Screw Dia25.4mm Microscope Accessories C - CS Converter A55.2201-B

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC
- Certification: CE
- Model Number: A55.2201-B
- Minimum Order Quantity: 1
- Packaging Details: Carton
- Supply Ability: 1000



Product Specification

- Product Name: Microscope Accessory C-CS Converter
- Inner Screw Dia.: 25.4mm
- Height: 5mm
- Model: A55.2201-B
- Theroy: Microscope Accessory
- Out Screw Dia.: 25.4mm
- Highlight: microscope accessory, parts of the microscope



Product Description

Inner And Out Screw Dia25.4mm Microscope Accessories C - CS Converter A55.2201-B

Specification

Model	A55.2201-B
Inner Screw Dia.	25.4mm
Out Screw Dia.	25.4mm
Height	5mm
Application	Microscope Accessory



A55.2201-B C-CS Converter



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China