



## Cool Light Fiber Optic Microscope Led Illuminator A56.2617



### Our Product Introduction

#### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A56.2617
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



#### Product Specification

- Model: A56.2617
- Power: 15W LED Lamp
- Product Name: Cool Light Fiber Optic LED Illuminator
- Packing Size: 285x230x255mm
- Size: 230x101x150mm
- LED Working Life: 50000H
- Highlight: cool light microscope led illuminator ,  
fiber optic microscope led illuminator ,  
led fiber optic microscope illuminator

## Product Description

### Cool Light Fiber Optic LED Illuminator Microscope Accessory A56.2617


- For observing temperature-sensitive specimens
- Gooseneck lights can be positioned independently
- Forced-air cooling system prevents overheating

#### Specification

A56.2617 Cool Light Fiber Optic LED Illuminator		
Power	15W LED Lamp	
Input voltage	10V-240V, 50Hz-60Hz	
CT	6000K	
Illuminance	Approximately 90000Lx	
Cooling System	Built In Large-area Radiator and Cooling Fan	
LED Working Life	50000H	
Fiber Optic Interfacetive	16mm (5/8")	
Size	230x101x150mm	
Net Weight	2.4kgs	
Gross Weight	2.9kgs	
Packing Size	285x230x255mm	
Optional Accessories		Item No.
Fiber Optic Light Guide (Gooseneck type)	Metal-sheathed Fiber Optic Bifurcated Light Guides with Condenser, Fiber :Length 550mm,Fiber Diameter 8mm, Standard Interface 5/8"	A56.2610-03
	Metal-sheathed Fiber Optic Single Light Guide with Condenser, Fiber Length 550mm,Fiber Diameter 8mm, Standard Interface 5/8"	A56.2610-04
	Metal-sheathed Fiber Optic Ring Light Guide with Condenser, Fiber Length 550mm,Fiber Diameter 8mm, Standard Interface 5/8"	A56.2610-0



A56.2617  
Cool Light Fiber Optic  
LED Illuminator







**Opto-Edu (Beijing) Co., Ltd.**



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China