



Magnification Digital Comparison Forensic Comparison Microscope A18.1808-C 2.7x~255x

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A18.1808-C
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: School|Laboratory|Research
- Adapter: SLR Camera|CCD
- Separation Line: Separation Line Adjusting System
- Illumination: LED
- Magnification: 2.7x~255x
- Theory: Comparison Microscope
- Highlight: **binocular optical microscope,**
bullet comparison microscope

Product Description

2.7x~255x Magnification Digital Comparison Microscope A18.1808-C

Specifications

A18.1808-C Digital Comparison Microscope

Magnification	2.7x~255x
Eyepiece	WF10X/22mm
	WF20X/14mm
Objective	Zoom Objective 0.7~4.9x, Zoom Ratio 1:7
	Auxiliary Lens 0.3x
	Auxiliary Lens 2.0x
Separation Line	Separation Line Adjusting System, Can Adjust Separation Line Easily
Magnification Correcting System	Magnification Correcting Lens 1 Pair, To Manually Adjust 2 Optical Systems Magnification To Be Exactly Same For Comparison
Bullet Holder	Manual Bullet Holder Set To Hold All Sizes Bullets, Shell Cartridge And Etc.
	Motorized Bullet Holder (3-speed)
Working Stage	Tilttable Big Working Stage, X-Y Manual Movement
	Motorized Focusing System, Z Axel Manual/Motorized Movement +Adjustable Arm For Holding Laptop
Light Source	High Intensity 150W Fiber Optic Cold Illumination System With Rigid Fiber Optic Cables. Includes Focused Collector Lens On Tips For Quick Adjustment Of Light Spot
	LED Multiband Cold Light Source With Blue, Green, Red, And White Colors. Colors Are Pushbutton Selected On Control Panel.
	5W LED Cold Light Illuminator With Manually Inserted Red And Green Color Filters
	Coaxial Light Lens
	Polarizing Lens
	Red/Green Filter
Adapter	SLR Camera Adapter
	CCD Adapter 0.5x, 1.0x
Control Box	Combined External Control Box
	Control Box Panel With Centralized Control Of Microscope Components. Controls Fiber Optic Light On/Off And Illumination Intensity, Controls LED Multiband Cold Light On/Off And Selection Of Four Colors, Controls On/Off Of LED Ring Lights, Controls 3-Speed Motorized Bullet Holder, Controls Up/Down Movement Of Motorized Stage Focusing
Software	Professional Forensic Image Analysis Software Crime Image
CCD	Built-in Digital CCD Built-in, 3.0M, Video/USB Output



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China