



## A23.0903-BL9 0.67x-4.5x Stereo Zoom Microscopes With Up & Bottom LED Illumination

Our Product Introduction

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A23.0903-BL9
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Applications: Industry
- Certification: CE|Rohs
- Eyepieces: WF10x/22mm, Both With Eye-cup
- Illumination: Up & Bottom LED Light
- Product Name: Stereo Zoom Microscope
- Magnification: 6.7~45x, Up To 2~225x With Auxiliary Lens
- Highlight: **stereo dissecting microscope,**  
**stereo binocular microscopes**

for more products please visit us on [cnoec.com](http://cnoec.com)

## Product Description

### A23.0903-BL9 0.67x-4.5x Zoom Stereo Microscope With Up & Bottom LED Illumination

#### Specifications:

A23.0903-BL9 Zoom Stereo Microscope, 0.67x-4.5x Specification		
Head	Brinocular, 45° Inclined, 360° Rotatable, +/-5° Diopter Adjustable, Interpupillary Adjustment 54~75mm	
Zoom Range	0.67~4.5x	
Zoom Ratio	1:6.5	
Eyepiece	WF10x/22mm, Both with Eye-cup	
Magnification	6.7~45x, Up to 2~225x with Auxiliary Lens	
Working Distance	100mm, With Standard 1x Lens. Up to 26~287mm With Auxiliary Lens	
Stand	Big Fan Track Stand BL9, With Focusing Knob, Focusing Range 50mm	
Illumination	Up & Bottom LED Light	
Power	100~240V Wide Range Voltage	
A23.0903-BL9 Zoom Stereo Microscope, 0.67x-4.5x Optional Accessories		
Eyepiece	WF15x/15mm	A51.0923-15
	WF20x/10mm	A51.0923-20
	WF25x/9mm	A51.0923-25
Auxiliary Lens	0.3x, W.D. 287mm	A52.0928-03
	0.4x, W.D. 217mm	A52.0928-04
	0.5x, W.D. 177mm	A52.0928-05
	0.75x, W.D. 117mm	A52.0928-75
	1.5x, W.D. 47mm	A52.0928-15
	2.0x, W.D. 26mm	A52.0928-20



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



[sale@optoedu.com](mailto:sale@optoedu.com)



[cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China