



WF10X Double Layers Mechanical Stage Biological Microscope A11.1316 With LED Lamp 1W Light Source

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1316
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Student|school|Laboratory|Hobby
- Certification: CE|Rohs
- Condenser: Abbe Condenser N.A.1.25
- Illumination: LED Lamp 1W
- Magnification: 1000X
- Product Name: Biological Microscope
- Highlight: **binocular biological microscope, biological compound microscope**



More Images



Product Description

WF10X Double Layers Mechanical Stage Biological Microscope A11.1316 With LED Lamp 1W Light Source

Specifications:

A11.1316 Biological Microscope		M	B	T
Head	Monocular Head, 30° Inclined, 360° Rotatable	•		
	Sliding Binocular Head, Interpupillary Distance:55-75mm, 45° Inclined, 360° Rotatable, Left & Right Diopter Adjustable		•	
	Sliding Trinocular Head, Interpupillary Distance:55-75mm, 45° Inclined, 360° Rotatable, Left & Right Diopter Adjustable			•
	Wide Field Plan WF10X			
Eyepiece	Wide Field Plan WF10X	•	•	•
	Achromatic 4X	•	•	•
	Achromatic 10X	•	•	•
	Achromatic 40X(S)	•	•	•
Objective	Achromatic 100X(S,Oil)	•	•	•
	Quadruple Nosepiece	•	•	•
	Stage	•	•	•
	Double Layers Mechanical Stage 140X155mm/75X50mm	•	•	•
Condenser	Abbe Condenser N.A.1.25, Φ2-Φ30mm Iris Diaphragm,With Φ32 Filter Holder	•	•	•
Focusing	Coaxial Coarse & Fine Focusing Adjustable Mechanism 25mm,Precision 0.002mm	•	•	•
Light Source	LED Lamp 1W, Brightness Adjustable	•	•	•
Weight	Gross Weight 7Kg, Net Weight 6Kg	•	•	•
Package	1 Set/ Carton, 33cm*28cm*43cm, Wood Cabinet Is Optional	•	•	•
A11.1316 Optional Accessory				
Eyepiece	Wide Field Plan WF16X			A51.1303-16



Head And Eyepiece Of A11.1316 Biological Microscope





Objective Of A11.1316 Biological Microscope



Package Of A11.1316 Biological Microscope



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China

