OPTO-EDU A59.5124 7"/10.1" 4K LCD Digital Camera 8.0M USB HDMI WIFI TF Card

Our Product Introduction


Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A59.5124
- Minimum Order Quantity: 1 pc
- Price: FOB $1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T,West Union,Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Resolution: Photo 5M/8M/12M Video 1080p
- Sensor: 4K 8.0M CMOS 1/2.8"
- Pixel Size: 1.45μm X 1.45μm
- Basic Function: Capture Photo, Record Video, Output To Computer +TF Card
- Highlight: **4K LCD digital microscope camera**, **USB HDMI microscope camera**, **WIFI TF card microscope camera**



for more products please visit us on cnoec.com

Product Description

7"/10.1" 4K LCD Screen With Built-in Measure Menu, 8.0M CMOS 1/2.8"
Output Type-C UVC+HDMI+2.4G WIFI+TF, PC Software StrangeView, Andorid/IOS APP T-Cam
Photo Resolution 5M/8M/12M, Video Resolution 1080p, Capture & Saved In TF Card
Including 18650 Rechargeable Battery (Not For Air Shipment) & 8G TF Card
Optional Adapter: -E = Eyepiece Adapter (Dia.23.2mm), -C = C-Mount Adapter (1" 32UN Screw)

OPTO-EDU (BEIJING) CO., LTD.
F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China
Tel:+8610 88696020 Fax:+8610 88696085

A59.5124

7"/10.1" 4K LCD Digital Camera, 8.0M, USB+HDMI+WIFI+TF Card



A59.5124 Specification

OPTO-EDU



7"/10.1" 4K LCD Screen 8.0M

Output Type-C UVC+HDMI+2.4G WIFI+TF+USB3.0

With Built-in Measure Menu

PC Software Strange View

Android / IOS APP T-Cam



	Button LCD Camera	2K Button LCD Camera	4K LCD Camera

for more products please visit us on cnoec.com

Model	A59.5121-			A59.5122-			A59.5124-	
	7E/7C	9E/9C	10E/10C	7E/7C	9E/9C	10E/10C	7E/7C	10E/10C
LCD Size	**7"**	**9"**	**10.1"**	**7"**	**9"**	**10.1"**	**7"**	**10.1"**
LCD Resolution	1024*600						1280*800	
Resolution	Photo 2M/5M/8M Video 1080p			Photo 2M/5M/8M/12M Video 1080p			Photo 5M/8M/12M Video 1080p	
Sensor	2.0M CMOS1/2.8"			2K 4.0M CMOS 1/3"			4K 8.0M CMOS 1/2.8"	
Pixel Size	2.9μm x 2.9μm			2μm x 2μm			1.45μm x 1.45μm	
Basic Function	**Capture Photo, Record Video, Output To Computer +TF Card**							
WIFI	-			2.4G, App T-Cam			2.4G/5G, App T-Cam	
HDMI	-			1080P-HDMI To Monitor			4K-HDMI To Monitor	
Mouse Control	-			Built-in Measure Menu			Built-in Adv. Measure Menu	
Type-C	Type-C Power In / USB-UVC Output To Computer, Measure Software Strange View							
USB	-			USBx1			USBx2	
TF Card	8G TF Card Support Max 64G						32G TF Card Support Max 64G	
Battery	18650-3.7V/2000MAH							
MINI USB	-			Connect AHD External Lens			-	
Adapter	**-E** = Eyepiece Adapter (Dia.23.2mm), **-C** = C-Mount Adapter (1" 32UN Screw)							
Power	Type-C Power In DC5V-2A							



Opto-Edu (Beijing) Co., Ltd.

0086 13911110627 sale@optoedu.com cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China